

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	12849-1KG-HON
<b>Description</b>	Copper(II) sulfate pentahydrate

---

## Technical specifications :

- ✓ copper(ii) sulfate 5h2o puriss meets bp, fcc, ph. eur. usp anal. spec. 1kg
- ✓ CAS N°7758-99-8
- ✓ Molecular weight 249.69 g/mol
- ✓ Molecular formula  $\text{CuO}_4\text{S} \cdot 5\text{H}_2\text{O}$
- ✓ Density 2.28 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form solid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1kg
  
- ✓ GHS pictograms GHS09 ; GHS07
- ✓ H Statements H302 ; H315 ; H319 ; H410
- ✓ P Statements P280

## Sales unit :

1kg Copper(II) sulfate pentahydrate

---