



Certificate of Compliance

No. 3N210713.BSPTU27

Test Report / Technical Construction File no. TCF-03-2020813-LVD+EMC

Certificate's
Holder:

Baoding Shenchen Precision Pump
Co., Ltd.
No.103, Building 2, Zhidian Industrial Park, Fuxing
East Road 999, Baoding, China

Certification ECM
Mark:



Product:
Model(s):

Peristaltic Pump
(see the following annex)

Verification to:

Standard:
EN 61010-1:2010+A1:2019, EN 61326-1:2013,
EN 61000-3-2:2019,
EN 61000-3-3:2013+A1:2019

related to CE Directive(s):
2014/35/EU (Low Voltage)
2014/30/EU (Electromagnetic Compatibility)

Remark: This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the technical documentation received from the manufacturer is satisfactory for the requirements of the ECM Certification Mark. The conformity mark above can be affixed on the products accordingly to the ECM regulation about its release and its use.

Additional information and clarification about the Marking:



The manufacturer is responsible for the CE Marking process, and if necessary, must refer to a Notified Body. This document has been issued on the basis of the regulation on ECM Voluntary Mark for the certification of products. RG01_ECM rev.3 available at: www.entecerma.it

Issuance date: 13 July 2021

Expiry date: 12 July 2026

Reviewer
Technical expert
Amanda Payne

Approver
ECM Service Director
Luca Bedonni

Ente Certificazione Macchine Srl

Via Ca' Bella, 243 – Loc. Castello di Serravalle – 40053 Valsamoggia (BO) - ITALY
☎ +39 051 6705141 📠 +39 051 6705156 ✉ info@entecerma.it 🌐 www.entecerma.it



Annex I

No. 3N210713.BSPTU27

Test Report / Technical Construction File no. TCF-03-2020813-LVD+EMC

Model(s):

LabS1, LabS3, LabK1, LabK3, LabQ, LabQ1, LabQ3, SK-Handypump, EM-UD15, EM-Handypump, EM-MiniPump, SP-MiniPump, Lab2015, BT100M, BT300M, BT600M, BT100N, BT300N, BT600N, BT100FJ, BT300FJ, BT600FJ, V1, V3, V6, V1-II, V3-II, F1, F1-II, F3, F3-II, F6, LabM1, LabM1-II, LabM1-III, LabM1-IV, LabM3, LabM3-II, LabM3-III, LabM3-IV, LabM6, LabM6-III, LabM6-IV, LabN1, LabN1-II, LabN1-III, LabN1-IV, LabN3, LabN3-II, LabN3-III, LabN3-IV, LabN6, LabN6-III, LabN6-IV, LabV1, LabV1-II, LabV1-III, LabV1-IV, LabV3, LabV3-II, LabV3-III, LabV3-IV, LabV6, LabV6-III, LabV6-IV, LabVU, LabF1, LabF1-II, LabF1-III, LabF1-IV, LabF3, LabF3-II, LabF3-III, LabF3-IV, LabF6, LabF6-III, LabF6-IV, LabGV, LabGM, LabGN, LabGF, LabPV-L, LabPV-M, LabPV-H, LabPF-L, LabPF-M, LabPF-H, LabPM-L, LabPM-M, LabPM-H, ZG600x, M6-3L, M6-6L, M6-12L, N6-3L, N6-6L, N6-12L, V6-3L, V6-6L, V6-12L, F6-3L, F6-6L, F6-12L, IF3, DF600, DF600Plus, DF600II, DF600PlusII, DF600III, DF600PlusIII, DF600IV, DF600PlusIV, DF350, DF350Plus, CF600, CF600Plus, CF600II, CF600PlusII, CF600III, CF600PlusIII, CF600IV, CF600PlusIV, CF350, CF350Plus, KF300, J25, EXP600, EXP600-A-I, EXP60-A-II, EXP600-B-I, EXP600-B-II, Exd600-3L, Exd600-6L, Exd600-12L, QD600, PDS, PDSII, PDSIII, PDSIV, BioD-I, BioD-II, OCMI, OCMI, OCMI, YZ1515x, YZ2515x, Minipump01, Minipump02, Handypump01, Handypump02, Micropump, MC series, AMC series, AMD series, AMDI, AMDII, AMDIII, AMDIV, AMDV, AMDVI, DG series, FD15, Deftpump, MT12, MT12II, MT12III, MD series, UC15, UC20, UC25, UD15, DY15, DY25, DZ25-3L, DZ25-6L, YZ35, SN25, SN15-16, SN15-14, DZ45, DZ45-I, DZ45-II, EasypumpI, EasypumpII, EasypumpIII, EasypumpIV, EasypumpV, EasypumpVI, EasypumpVII, EasypumpVIII, EasypumpIX, EasypumpX