

TECHNICAL DATASHEET**PART NUMBER****130340CL**

The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.

**ClearLine Tips**

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

| | |
|----------------------|------------------------------|
| Part number | 130340CL |
| Packaging | ClearLine rack with slip lid |
| Volume | 10 µl |
| Graduation | 1, 2, 5, 10 |
| Length | 44,8 mm |
| Color | Natural |
| Filter | Yes |
| Low Retention | No |
| Sterility | Yes |
| Autoclavable | No |

Sales units :

960 tips