



ClearLine® Culture Flask

Description : *Sterile tissue culture flasks.*

Features :

- ✓ Flasks made in crystal polystyrene, cap in polyethylene
- ✓ Packaged in zip resalable bags
- ✓ Manufactured in ISO 8 Cleanroom conditions
- ✓ Gamma sterilized
- ✓ RNase, DNase-free, endotoxins certified <0.5EU/ml
- ✓ Lot n° and expiry date on each flask, certificates available upon request
- ✓ Graduations are engraved and visible
- ✓ Easy stacking
- ✓ Writing area

Volume (ml)	Cell growth surface (cm ²)	Units /bag	Units /box	Cap
50	25	10	200	filtered