

**TECHNICAL DATASHEET**

**REFERENCE :**



**ClearLine®  
Centrifuge Tubes**

**Description :**

*Polypropylene centrifuge tubes,  
Available sterile or non-sterile,  
15 ml & 50 ml  
Racked or in zip resealable*

**Features :**

- ✓ Tubes in virgin medical grade polypropylene (PP)
- ✓ Caps in high-density polyethylene (HDPE)
- ✓ Excellent transparency
- ✓ Max RCF 12000 x g
- ✓ Autoclavable at 121°C (20 min), can be stored down to -80°C
- ✓ Frosted writing area
- ✓ Graduated
- ✓ Manufactured in Cleanroom ISO Class 8 conditions
- ✓ **Gamma sterilised**
- ✓ **RNase, DNase-free, endotoxins certified < 0,5 EU /ml**
- ✓ **Lot n° & expiration date on each box**, certificates available upon request

Cat No	Volume (ml)	Sterile	Packaging	Units /box
131032C	50	Yes	Bags of 25	500
131135C	50	Yes	Racks of 25	500
131034C	50	Non	Bags of 25	500
131036C	50	Yes	Bags of 25 skirted tubes	500
131038C	15	Yes	Racks of 25	500
131040C	15	Non	Bags of 25	500
131042C	15	Yes	Bags of 25	500