

**TECHNICAL DATASHEET**

**REFERENCE :**



**ClearLine®  
Centrifuge Tubes**

**Description :**

*Polypropylene centrifuge tubes,  
Available sterile or non-sterile,  
15 ml & 50 ml  
Racked or in zip resealable*

**Features :**

- ✓ Tubes in virgin medical grade polypropylene (PP)
- ✓ Caps in high-density polyethylene (HDPE)
- ✓ Excellent transparency
- ✓ Max RCF 12000 x g
- ✓ Autoclavable at 121°C (20 min), can be stored down to -80°C
- ✓ Frosted writing area
- ✓ Graduated
- ✓ Manufactured in Cleanroom ISO Class 8 conditions
- ✓ **Gamma sterilised**
- ✓ **RNase, DNase-free, endotoxins certified < 0,5 EU /ml**
- ✓ **Lot n° & expiration date on each box**, certificates available upon request

| Cat No  | Volume (ml) | Sterile | Packaging                | Units /box |
|---------|-------------|---------|--------------------------|------------|
| 131032C | 50          | Yes     | Bags of 25               | 500        |
| 131135C | 50          | Yes     | Racks of 25              | 500        |
| 131034C | 50          | Non     | Bags of 25               | 500        |
| 131036C | 50          | Yes     | Bags of 25 skirted tubes | 500        |
| 131038C | 15          | Yes     | Racks of 25              | 500        |
| 131040C | 15          | Non     | Bags of 25               | 500        |
| 131042C | 15          | Yes     | Bags of 25               | 500        |