



**ClearLine® Tissue Culture Dish**

**Description :** Tissue Culture treated polystyrene sterile Petri dish

**Features :**

- ✓ Tissue culture treated (adherent)
- ✓ Ultra-flat bottom and uniform thickness across
- ✓ Culture surface is smooth and unstriated
- ✓ 4 venting slots ensure optimal gaseous exchange
- ✓ Packaged in zip resealable self-standing plastic bags
- ✓ Lot n° on each bag, certificates are available on demand
- ✓ **Sterilised by gamma irradiation**
- ✓ **Certified RNase, DNase-free and endotoxins are < 0,5 EU /ml**

Cat No	Petri dish type	Diameter* (mm)	Dimensions** (h x Ø mm)	Cell-growth area (cm <sup>2</sup> )	Units /bag	Units /box
131050C	100 mm	87,8	22,0 x 95,6	60,8	10	300

\* inner diameter ; \*\* with closed lid