

TECHNICAL DATA SHEET

Part Number	13440-500ML-HON
Description	Sodium hypochlorite solution

Technical specifications :

- ✓ sodium hypochlorite solution, 6-14% active chlorine basis 500ml
- ✓ CAS N°7681-52-9
- ✓ Chemical purity 6-14%
- ✓ Molecular weight 74.44 g/mol
- ✓ Molecular formula ClNaO
- ✓ Density 1.206 g/cm³ (25 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition 2-8°C
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL

- ✓ GHS pictograms GHS09 ; GHS05
- ✓ H Statements H410 ; H314 ; H290
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

500mL Sodium hypochlorite solution
