

## TECHNICAL DATASHEET

## PART NUMBER

## 134750CL




The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



### ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

|               |                              |
|---------------|------------------------------|
| Part number   | 134750CL                     |
| Packaging     | ClearLine rack with slip lid |
| Volume        | 1250 µl                      |
| Graduation    | 200, 500, 1000               |
| Length        | 102 mm                       |
| Color         | Natural                      |
| Filter        | No                           |
| Low Retention | No                           |
| Sterility     | No                           |
| Autoclavable  | Yes                          |



### Sales units :

768 tips