

## Description

- Reagent bottle Wide-Mouth

## Part number

- 135745

## Pieces by sales :

1000 Bottles



---

## Description

Suitable for storage, transportation and packing of liquid, with uniform thickness of wall, sturdy and durable

## Technical data

- Produced in High Density Polyethylene (HDPE) with liner less polypropylene caps to ensure leakproof performance and loose material contamination
- Wide top design for easy filling and reduction of eddy currents caused by perfusion
- Chemical resistance to most acid, alkali and alcohol
- Volume: 15 ml
- Color: brown