

TECHNICAL DATA SHEET

Part Number	13767-1KG-HON
Description	Silica gel orange

Technical specifications :

- ✓ silica gel orange, with moisture indicator free of heavy metals 1kg
- ✓ CAS N°112926-00-8
- ✓ Molecular formula SiO₂
- ✓ Density 2 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form Beads

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1kg

- ✓ P Statements P280

Sales unit :

1kg Silica gel orange
