

Guest Medical

Unit A6, Larkfield Trading Estate,
Aylesford, Kent, ME20 6SW, United Kingdom
Telephone : +44 (0) 1622 791895
email: enquiries@guest-medical.co.uk
www.guest-medical.co.uk

B8617 Medium Autoclave Bags

Product specification

- Material:** High temperature clear polypropylene polymer
- Size:** 420mm x 600mm
- Thickness:** 38 microns
- Print Artwork:** Biohazard symbol, Hazardous Waste text
- Print Artwork Colour:** Pantone 286
- Temperature for use:** 134°C - 145°C (maximum temperature)
- Quantity per carton:** 500 bags
- Carton Size:** 483mm x 340mm x 102mm
- Carton weight:** 9.35Kg

Signed. 

Dated: 14th December 2020

Natalie Ferguson – Microbiologist & Technical Support Manager
For and on behalf of Guest Medical Limited