

FICHE TECHNIQUE

Référence	Voir tableau
Désignation	Ruban adhésif Time tape



Description :

- ✓ Longueur 12 m
- ✓ Diamètre intérieur 26 mm

Largeur (mm)	12	19	25	38	50
Blanc	140259A	140250A	140251A	140252A	140253A
Bleu	140259B	140250B	140251B	140252B	140253B
Cuivre	140259C	140250C	140251C	140252C	140253C
Rose	140259D	140250D	140251D	140252D	140253D
Tilleul	140259E	140250E	140251E	140252E	140253E
Fuschia	140259F	140250F	140251F	140252F	140253F
Gris	140259G	140250G	140251G	140252G	140253G
Chartreuse	140259H	140250H	140251H	140252H	140253H
Jaune	140259J	140250J	140251J	140252J	140253J
Lavande	140259L	140250L	140251L	140252L	140253L
Violet	140259M	140250M	140251M	140252M	140253M
Noir	140259N	140250N	140251N	140252N	140253N
Orange	140259P	140250P	140251P	140252P	140253P
Ocre	140259Q	140250Q	140251Q	140252Q	140253Q
Rouge	140259R	140250R	140251R	140252R	140253R
Saumon	140259S	140250S	140251S	140252S	140253S
Turquoise	140259T	140250T	140251T	140252T	140253T
Vert	140259V	140250V	140251V	140252V	140253V

Unité de vente :

1