

TECHNICAL DATA SHEET

| | |
|--------------------|--------|
| Part Number | 140262 |
| Description | Tape |



Benefits :

- ✓ Set of 12 rolls of Time Tape colored adhesive tape
- ✓ Assorted colors : 2 white, 2 yellow, 2 green, 2 blue, 2 red, 1 pink and 1 orange
- ✓ Autoclavable (+130°C) and usable in the freezer (-50°C)
- ✓ Can be used in ovens
- ✓ Stick and peel without leaving a trace
- ✓ Marking in pencil, ink or felt-tip pen

Technical specifications :

- ✓ Length of ribbons : 12 m
- ✓ Width of ribbons : 25 mm

Sales unit :

1 Set of 12 rolls of colored adhesive tape