



sharpsafe
SHARPSAFE 5^{ÈME} GENERATION 

General data

Item code	SHA-51652610
GTIN code	5056271508320
Product Name	Sharpsafe® 5 th GEN. 0,8 liter
Description	Sharps container 0,8 liter composed by 1 body + 1 lid to assemble
Brand	Sharpsafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste

Characteristics data

Nominal volume	0,8 l
Total volume	0,79 l
Usefull volume	0,67 l
Maximum weight	0,31 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
Dimensions (assembled) in mm :	
Long	108 mm
Wide	90 mm
Height	149 mm
Opening	58 x 36 mm
Fill line height	122 mm
Weight	103 g
Material	Polypropylene recycled 100% (body) + Polypropylene virgin 80% & recycled 20% (lid)
Phtalates & Latex free	Yes
Carbon footprint reduction	73%
Autoclavable	Yes
Autoclavability conditions	134° during 18 mn with lid open
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availability Accessories	Stands, trays, trolley brackets (see brochure and compatibility tables)


Logistic data

Carton quantity	120
Carton weight	12,7 Kg
Carton dimensions (LxWxH)	540 x 380 x 428 mm
Carton volume	87,8256 dm ³
Palet quantity	2400 units
Pallet weight	269 Kg
Pallet dimension (LxWxH)	120 x 80 x 230 cm

Certifications data

ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark - (NF 302 in progress)
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes

Manufacture data

Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group