



**sharpsafe**  
SHARPSAFE  
5<sup>ÈME</sup> GENERATION 

General data	
Item code	SHA-51602610
GTIN code	5056271508337
Product Name	Sharpsafe® 5 <sup>th</sup> GEN. 1 liter
Description	Sharps container 1 liter componed by 1 body + 1 lid to assemble
Brand	Sharpsafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste
Characteristics data	
Nominal volume	1 l
Total volume	1,03 l
Usefull volume	0,85 l
Maximum weight	0,31 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
Dimensions (assembled) in mm :	
Long	108 mm
Wide	90 mm
Height	187 mm
Opening	58 x 36 mm
Fill line height	151 mm
Weight	123 g
Material	Polypropylene recycled 100% (body) + Polypropylene virgin 80% & recycled 20% (lid)
Phtalates & Latex free	Yes
Carbon footprint reduction	64%
Autoclavable	Yes
Autoclavability conditions	134° during 18 mn with lid open
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availability Accessories	Stands, trays, trolley brackets (see brochure and compatibility tables)
Logistic data	
Carton quantity	100
Carton weight	13,7 Kg
Carton dimensions (LxWxH)	540 x 380 x 428 mm
Carton volume	87,8256 dm <sup>3</sup>
Palet quantity	2000 units
Pallet weight	289 Kg
Pallet dimension (LxWxH)	120 x 80 x 230 cm
Certifications data	
ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark - (NF 302 in progress)
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes
Manufacture data	
Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group

