



**sharpsafe**  
SHARPSAFE 5<sup>ÈME</sup> GENERATION 

**General data**

Item code	SHA-50502610
GTIN code	5056271508344
Product Name	Sharpsafe® 5 <sup>th</sup> GEN. 1,8 liter
Description	Sharps container 1,8 liter composed by 1 body + 1 lid to assemble
Brand	Sharpsafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste

**Characteristics data**

Nominal volume	1,8 l
Total volume	1,75 l
Usefull volume	1,45 l
Maximum weight	0,4 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
Dimensions (assembled) in mm :	
Long	124 mm
Wide	114 mm
Height	220 mm
Opening	57 x 53 mm
Fill line height	178 mm
Weight	169 g
Material	Polypropylene recycled 100% (body) + Polypropylene virgin 80% & recycled 20% (lid)
Phtalates & Latex free	Yes
Carbon footprint reduction	67%
Autoclavable	Yes
Autoclavability conditions	134° during 18 mn with lid open
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availability Accessories	Stands, trays, trolley brackets (see brochure and compatibility tables)


**Logistic data**

Carton quantity	75
Carton weight	12,8 Kg
Carton dimensions (LxWxH)	606 x 364 x 414 mm
Carton volume	91,3218 dm <sup>3</sup>
Palet quantity	1500 units
Pallet weight	271 Kg
Pallet dimension (LxWxH)	120 x 80 x 230 cm

**Certifications data**

ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark - (NF 302 in progress)
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes

**Manufacture data**

Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group