



sharpsafe
SHARPSAFE
5^{ÈME} GENERATION 

General data	
Item code	SHA-51406610 (without flap) - SHA-51405610 (with flap)
GTIN code	5056271508368 (without flap) - 5056271508351 (with flap)
Product Name	Sharpsafe® 5 th GEN. 2 liters (flap in option)
Description	Sharps container 2 liters composed by 1 body + 1 lid to assemble
Brand	Sharpsafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste
Characteristics data	
Nominal volume	2 l
Total volume	2,23 l
Usefull volume	1,67 l
Maximum weight	0,5 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
Dimensions (assembled) in mm :	
Long	197 mm
Wide	126 mm
Height	177 mm
Opening	55 x 55 mm
Fill line height	128 mm
Weight	220 g (228 g with flap)
Material	Polypropylene recycled 100% (body) + Polypropylene virgin 80% & recycled 20% (lid)
Phtalates & Latex free	Yes
Carbon footprint reduction	61%
Autoclavable	Yes
Autoclavability conditions	134° during 18 mn with lid open
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availability Accessories	Stands, trays, trolley brackets (see brochure and compatibility tables)
Logistic data	
Carton quantity	50
Carton weight	12,6 Kg
Carton dimensions (LxWxH)	585 x 385 x 375 mm
Carton volume	84,459375 dm ³
Palet quantity	1200 units
Pallet weight	318 Kg
Pallet dimension (LxWxH)	120 x 80 x 230 cm
Certifications data	
ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark - (NF 302 in progress)
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes
Manufacture data	
Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group

