


BagPage® + Patented Full-page filter bag for blender

Ref 121 025 (100 mL) - Ref 122 025 (400 mL) - Ref 122 200 (2000 mL) - Ref 123 010 (3500 mL)



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagPage® + is a full-page filter blender bag with a full surface micro-perforated filter. Bacteria are extracted from the sample and goes through the filter during blending. You can then pipet behind the filter with no residue or particles of the sample.

TECHNICAL SPECIFICATIONS

- Bag with full-page filter
- Ideal for homogenizing pasty samples
- **Multilayer®**: reinforced multicoated complex
- Full-page calibrated microperforated filter
- Filter porosity: 280 microns
- Rigid and transparent
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenization of the sample
- Compatible with any lab blender
- Approved for food contact: Regulation (EC) No 1935/2004
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Average shelf-life: 25 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 100/400/2000/3500 mL
- In compliance with ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

BagPage® 100

Ref 121 025

- Max blending volume: 100 mL
- Optimal blending volume: 5 - 50 mL
- Bag dimensions: 95 x 180 mm
- Box dimensions: 34 x 22 x 12 cm, weight: 2 kg
- Pack of 25
- Box of 500

BagPage® + 400

Ref 122 025

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 190 x 300 mm
- Box dimensions: 34 x 22 x 26 cm, weight: 6 kg
- Pack of 25
- Box of 500

BagPage® + 2000

Ref 122 200

- Max blending volume: 2000 mL
- Optimal blending volume: 400 - 1500 mL
- Bag dimensions: 250 x 380 mm
- Box dimensions: 47.5 x 28.5 x 13.5 cm, weight: 4.9 kg
- Pack of 25
- Box of 250

BagPage® + 3500

Ref 123 010

- Max blending volume: 3750 mL
- Optimal blending volume: 400 - 3750 mL
- Bag dimensions: 380 x 600 mm
- Box dimensions: 41 x 33.5 x 11 cm, weight: 4.2 kg
- Pack of 10
- Box of 100

BagPage® + 400 unmarked

Ref 122 000

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 190 x 300 mm
- Box dimensions: 34 x 22 x 26 cm, weight: 6 kg
- Pack of 25
- Box of 500