

# Technical Datasheet

---

Article No. 1411

Methanol technical grade  
min. 99.5 %

For Laboratory use.

| Parameter        | Value             |
|------------------|-------------------|
| Appearance       | colourless liquid |
| pH- value        | ± 7               |
| Melting Point    | - 98 °C           |
| Boiling Point    | 64.7 °C           |
| Flash Point      | 11 °C             |
| Density (20 °C)  | 0.79 g/ml         |
| Mol.-Weight      | 32.04 g/mol       |
| H <sub>2</sub> O | max. 0.15 %       |
| Free acids       | max. 0.003 %      |