

[Argos Technologies](#) – Item #04396-40



Specifications & Description

- Handle length (in) - 6
- Overall Length (in) - 6
- Handle Type - Brass
- Description - Brass Inoculating Loop and Needle Holder, 6"; Each

Reusable Inoculating Loops with handles are used to select, inoculate or transfer microbiological samples. Handles are insulated to protect the hand during use. Loops and Handles cool rapidly after sterilization. Loops and without handles are sold separately as models 04396-35 through -39. Needle without handle is available as model 04396-49.