



Declaration of conformity



According to the in Vitro Diagnostic Medicine Device Directive 98/79/EC

Manufacturer: GUANGZHOU JET BIO-FILTRATION CO., LTD
No.1 Doutang Road, Yonghe Economical & Technological, Development
District, Luogang, 511356 Guangzhou, Guangdong, P.R. China
No.173-175 Youyi Road Economical & Technological Development
District 510730 Guangzhou, P.R. China

Authorised

Representative: TÜV Süddeutschland Group
Ridler str.65, Munchen, 80339 Deutschland

Description/Category: In Vitro Diagnostic Medical Devices

IVDD Classification : Other Type (General)

Products Name : See attached list

The undersigned hereby declares that the in Vitro Diagnostic Medical Devices as specified in the attached list conform with the essential requirements listed in the Annex I of European In Vitro Diagnostic Medicine Device Directive 98/79/EC(IVDD).

This declaration of conformity is based on the European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD). Annex III is supported by a EC-REP TUV SUD Product Service GmbH Annex V Certification. With reference to articles 1 and 3 of the IVDD.

Quality Manager

Guangzhou
Aug. 17, 2015
(Place and date of issue)

York Zhang
(Name and signature of authorized person)



Declaration of Conformity Directive 98/79/EC

PRODUCTS LIST

| Product Type | Cat.No. | Description and Packing |
|---------------------------|--|--|
| Syringe Driven Filters | 23966 | 13mm, 0.22um, PES, Green, Non-sterile, 100/pk(box), 1000/cs |
| | 23967 | 13mm, 0.45um, PES, Green, Non-sterile, 100/pk(box), 1000/cs |
| | 23968 | 25mm, 0.22um, PES, Green, Non-sterile, 100/pk(box), 1000/cs |
| | 23969 | 25mm, 0.45um, PES, Green, Non-sterile, 100/pk(box), 1000/cs |
| | 23970 | 4mm, 0.22um, Nylon, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23971 | 4mm, 0.45um, Nylon, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23972 | 4mm, 0.22um, CA, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23973 | 4mm, 0.22um, CA, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23974 | 4mm, 0.22um, PTFE, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23975 | 4mm, 0.45um, PTFE, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23976 | 4mm, 0.22um, PVDF, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23977 | 4mm, 0.45um, PVDF, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23978 | 4mm, 0.22um, PES, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23979 | 4mm, 0.45um, PES, Clear, Non-sterile, 100/pk(box), 1000/cs |
| | 23980 | 30mm, 0.22um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs |
| | 23981 | 30mm, 0.45um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs |
| | 23982 | 30mm, 0.22um, Cellulose, Yellow, Non-sterile, 100/pk(box), 1000/cs |
| | 23983 | 30mm, 0.45um, CA, Red, Non-sterile, 100/pk(box), 1000/cs |
| | 23984 | 30mm, 0.22um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs |
| | 23985 | 30mm, 0.45um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs |
| | 23986 | 30mm, 0.22um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs |
| | 23987 | 30mm, 0.45um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs |
| | 23988 | 30mm, 0.22um, PES, Green, Non-sterile, 100/pk(box), 1000/cs |
| | 23989 | 30mm, 0.45um, PES, Green, Non-sterile, 100/pk(box), 1000/cs |
| | 26142 | 13mm, 0.22um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs |
| | 26143 | 13mm, 0.45um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs |
| 26144 | 25mm, 0.22um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs | |

| | | |
|-------------------------|---|--|
| 26145 | 25mm, 0.45um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs | |
| 26146 | 13mm, 0.22um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs | |
| 26147 | 13mm, 0.45um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs | |
| 26148 | 25mm, 0.22um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs | |
| 26149 | 25mm, 0.45um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs | |
| 26150 | 13mm, 0.22um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs | |
| 26151 | 13mm, 0.45um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs | |
| 26152 | 25mm, 0.22um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs | |
| 26153 | 25mm, 0.45um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs | |
| 26155 | 13mm, 0.45um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs | |
| 26156 | 13mm, 0.22um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs | |
| 26157 | 25mm, 0.45um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs | |
| 26158 | 25mm, 0.22um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs | |
| Serological Pipettes | 56400-11 | Pipet,Bact,Plas,Milk,1.1mL,Multi Pack,25/Pack,1000/Case |
| | 56400-22 | Pipette, Sero, Plas,TD, Multi. 2.2ml, 2x 1/100mL, 25/pack, 1000/cs |
| | 56404-11 | Pipet,Bact,Plas,Milk,1.1mL,Individually Wrapped, 200/Pack,1000/Case |
| | 56404-22 | Pipette, Sero, Plas,TD, 2.2ml, 2x 1/100mL, Individually Wrapped, 500/CS |
| | 56700-10110 | Pipet,Sero,Plas,TD,Canister Pack,10ml,Plugged,Sterile,50/Canister,500/Case |
| | 56700-11100 | Pipet,Sero,Plas,TD,Canister Pack,1ml,Plugged,Sterile,50/Canister,500/Case |
| | 56700-21100 | Pipet,Sero,Plas,TD,Canister,2ml,Plugged,50/Canister,500/Case |
| | 56700-5110 | Pipet,Sero,Plas,TD,Canister,5ml,Plugged,50/Canister,500/Case |
| | 56800-10110 | Pipette,Sero,Plas,TD,Multi.10mL,10x 1/10mL,Shelfpacks of 25,Case of 500 |
| | 56800-1110 | Pipette,Sero,Plas,TD,Multi.1mL,1x 1/10mL,Shelfpacks of 25,Case of 1000 |
| | 56800-11100 | Pipet, Sero, Plas, TD, Multi, 1ml, Plugged, Sterile, 25/Pack, 1000/Case |
| | 56800-21100 | Pipet,Sero,Plas,TD,Multi Pack,2ml, 25/Pack,500/Case |
| | 56800-25210 | Pipet,Sero,Plas,TD,Multi Pack,25mL,Plugged, Sterile,10/pack,200/case. |
| | 56800-50510 | Pipette,Sero,Plas,TD,Multi.50mL,50x 5/10mL,Shelfpacks of 10,Case of 100 |
| | 56800-5110 | Pipet,Sero,Plas,TD,Multi Pack,5ml,Plastic, Plugged, Sterile, 25/Pack,500/Case |
| | 56825-10110 | Pipet, Sero, Plas, TD, Multi Pack, 10mL,Plugged, Sterile, 25/Pack, 500/Case *Wide Tip* |
| | 56900-10010 | 100ml,Sterilized,Pink,Ind Wrap(Paper/Plastic),50/cs |
| | 56900-10110 | Pipette,Sero,Plas,TD,10ml,Individually Wrapped, Plugged, Sterile,200/Case (Shelfpacks of 50) |

| | | |
|---------------------|-------------|--|
| | 56900-1110 | Pipet,Sero,Plas,TD,1mL,Individually Wrapped, Plugged, Sterile,1000/Case (Shelfpacks of200) |
| | 56900-11100 | Pipet,Sero,Plas,TD,1ml,Individually Wrapped, Plugged, Sterile,1000/Case (Shelfpacks of 200) |
| | 56900-21100 | Pipet,Sero,Plas,TD,2ml,Individually Wrapped, Plugged, Sterile,500/Case (Shelfpacks of 125) |
| | 56900-25210 | Pipette,Sero,Plas,TD,25mL,Individually Wrapped, Plugged, Sterile,200/Case (Shelfpacks of 50) |
| | 56900-50510 | Pipet, Sero, Plas, TD, 50ml,Individually Wrapped, Plugged,Sterile, 100/Case (Shelfpacks of 25) |
| | 56900-5110 | Pipet,Sero,Plas,TD,5mL,Individually Wrapped, Plugged, Sterile,200/Case (Shelfpacks of 50) |
| Aspirating Pipets | APS10 | 10ml,Sterilized (Individually wrapped ,200/cs) |
| | APS2 | 2ml,Sterilized (Individually wrapped,500/cs) |
| | APS5 | 5ml,Sterilized (Individually wrapped,200/cs) |
| Reservoir Bottle | B0001 | 1000ml, Non-sterilized (1/pkg, 24/cs) |
| Pipetting Reservoir | B3110-100 | 100mL, Sterile,DNA-RNase & pyrogen free, 5/pk, 100pcs/case. |
| | B3110-100NS | 100mL, Bulk Pack , Non Sterile, DNA-RNase & pyrogen free, 100pcs/bag, 1bag/case. |
| | B3110-50 | 100mL, Sterile, Individually Wrapped, DNA-RNase & pyrogen free, 50pcs/case. |
| | B3125-100 | 25mL,Sterile,DNA-RNase & pyrogen free,5/pk, 100/case. |
| | B3125-100NS | 25mL , Non Sterile, DNA-RNase & pyrogen free, Bulk Pack, 100pcs/bag, 1bag/case. |
| | B3125-50 | 25mL, Sterile, DNA-RNase & pyrogen free, Individually Wrapped, 50pcs/case. |
| | B3150-100 | 50mL, Sterile,DNA-RNase & pyrogen free, 5/pk, 100pcs/case. |
| | B3150-100NS | 50mL, Bulk Pack , Non Sterile, DNA-RNase & pyrogen free, 100pcs/bag, 1bag/case. |
| | B3150-50 | 50mL, Sterile, DNA-RNase & pyrogen free, Individually Wrapped, 50pcs/case. |
| 12 Channel basin | B3134 | 12 Channel non-sterile basin, individually packed, 50/cs |
| | B3135 | 12 Channel sterile basin, Individually packed, 20/box, 12boxes/case |
| Vacuum Filtration | BMC2225 | 250ml, MCE, 0.22um, Sterilized (1/pkg, 12/cs) |
| | BMC2250 | 500ml, MCE, 0.22um, Sterilized (1/pkg, 12/cs) |
| | BMC4525 | 250ml, MCE, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BMC4550 | 500ml, MCE, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BNY2225 | 250ml,NYLON,0.22um, Sterilized (1/pkg,12/cs) |
| | BNY2250 | 500ml,NYLON,0.22um, Sterilized (1/pkg,12/cs) |
| | BNY4525 | 250ml, NYLON, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BNY4550 | 500ml, NYLON, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BPE2210 | 1000ml, PES, 0.22um, Sterilized (1/pkg, 12/cs) |
| | BPE2215 | 150ml, PES, 0.22um, Sterilized (1/pkg, 12/cs). |

| | | |
|---------------------|--------------|---|
| | BPE2225 | 250ml, PES, 0.22um, Sterilized (1/pkg,12/cs) |
| | BPE2250 | 500ml, PES, 0.22um, Sterilized (1/pkg, 12/cs) |
| | BPE4510 | 1000ml, PES, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BPE4515 | 150ml, PES, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BPE4525 | 250ml, PES, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BPE4550 | 500ml, PES, 0.45um, Sterilized (1/pkg, 12/cs) |
| Vacuum Filtration | BPV2210 | 1000ml, PVDF, 0.22um, Sterilized (1/pkg, 12/cs) |
| | BPV2225 | 250ml, PVDF, 0.22um, Sterilized (1/pkg, 12/cs) |
| | BPV2250 | 500ml, PVDF, 0.22um, Sterilized (1/pkg, 12/cs) |
| | BPV4510 | 1000ml, PVDF, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BPV4525 | 250ml, PVDF, 0.45um, Sterilized (1/pkg, 12/cs) |
| | BPV4550 | 500ml, PVDF, 0.45um, Sterilized (1/pkg, 12/cs) |
| Cell Scraper | C2500 | 25cm ,Cell Scraper,Purple,Sterilized,Non-pyrogenic,Lifter Blade(Individually wrapped, 100pcs per bag/case) |
| | C3900 | 39cm ,Cell Scraper,Sterilized,Non-pyrogenic,Lifter Blade(Individually wrapped, 100pcs per |
| Cell Strainer | CS040 | 40µm Cell Strainer, Blue, Sterile, Individually Packaged, 50/Case |
| | CS070 | 70µm Cell Strainer, White, Sterile, Individually Packaged, 50/Case |
| | CS100 | 100µm Cell Strainer, Yellow, Sterile, Individually Packaged, 50/Case |
| Pipetting Reservoir | D3110-100 | Denville Basins, 100mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free (100/case). Denville # 1005110. |
| | D3110-100NS | Denville Basins, 100mL, Bulk Pack, Non Sterile, 100/cs. Denville # 1005118 |
| | D3110-50 | Denville Basins, 100mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/Case. Denville#1005108. |
| | D3125-100 | Denville Basins, 25mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free 100/case. Denville # 1005120 |
| | D3125-100NS | Denville Basins, 25mL, Bulk Pack, Non-Sterile, 100/cs. Denville # 1005121 |
| | D3125-50 | Denville Basins, 25mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/case. Denville # 1005119 |
| | D3150-100 | Denville Basins, 50mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free 100/case. Denville # 1005109 |
| | D3150-100NS | Denville Basins, 50mL, Bulk Pack, Non Sterile, 100/cs. Denville # 1005111 |
| | D3150-50 | Denville Basins, 50mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/case. Denville # 1005122 |
| | DB3110-50 | ClearLine Basins, 100mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free 50/Case |
| | DB3125-100 | ClearLine Basins, 25mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free 100/case |
| | DB3125-100NS | ClearLine Basins, 25mL, Bulk Pack, Non-Sterile, 100/cs |
| | DB3125-50 | ClearLine Basins, 25mL, Sterile, Individually, Wrapped, DNA-RNase & Pyrogen Free 50/case |

| | | |
|----------------------|------------|--|
| | DB3150-50 | ClearLine Basins, 50mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/case |
| Cell Scraper | DC2500 | 25cm ,Cell Scraper, Purple, Sterilized, Nonpyrogenic, Lifter Blade (Individually wrapped, 100pcs per bag/case) |
| | DC3900 | 39cm ,Cell Scraper, Purple, Sterilized, Nonpyrogenic, Lifter Blade (Individually wrapped, 100pcs per case) |
| Tissue Culture Dish | DDT1000 | Tissue culture dishes, 100mm, surface treated, sterile, 10/bag, 300/case |
| | DDT3500 | Tissue culture dishes, 35mm, surface treated, sterile, 10/bag, 960/case |
| | DDT7000 | Tissue culture dishes, 70mm, surface treated, sterile, 10/bag, 600/case |
| Tissue Culture Flask | DF1105 | Tissue culture flask, 50ml, 25cm ² , non-vented, sterile with treated surface,10/bag 200/case |
| | DF1125 | Tissue culture flask, 250ml, 75cm ² ,non-vented, sterile with treated surface,5/bag ,100/case |
| | DF1160 | Tissue culture flask, 600ml, 182cm ² ,non-vented, sterile with treated surface,5/bag, 40/case |
| | DF1205 | Tissue Culture Flask, 50ml with Vented Cap, Sterile with Treated Surface,10/Bag, 200/Case |
| | DF1225 | Tissue Culture Flask, 250ml, with Vented Cap, Sterilize with Treated Surface, 5/Bag,100/Case. |
| | DF1260 | Tissue Culture Flask, 600ml with Vented Cap, Sterile with Treated Surface, 5/Bag 40/Case |
| | F1205 | 50ml, Surface-treated, w/vent cap, Sterilized (10/pkg, 200/cs) |
| Tissue Culture Plate | DP012096 | 96 well, U-bottom, Surface-treated, Sterilized, 1/pkg, 100/cases |
| | DP1006 | 6well, Surface-treated, Sterilized (1/pkg, 100/cs) |
| | DP1012 | 12 well, Surface-treated, Sterilized, 1/pkg, 100/cs |
| | DP1024 | 24 well, Surface-treated, Sterilized, 1/pkg, 100/case |
| | DP1096 | 96 well, Surface-treated, Sterilized, 1/pkg, 100/case |
| | P0096 | 96well, Non-treated, Sterilized (1/pkg, 100/cs) |
| | P1006 | 6well, Surface-treated, Sterilized (1/pkg, 100/cs) |
| | P1012 | 12well, Surface-treated, Sterilized (1/pkg, 100/cs) |
| | P1024 | 24well, Surface-treated, Sterilized (1/pkg, 100/cs) |
| | P1096 | 96well, Surface-treated, Sterilized (1/pkg, 100/cs) |
| | P2096 | 96well,without lid, non-treated and non-steriled,(10/bag,100/cs). |
| Centrifuge Tubes | 4503 | 50ml,Conical,Sterilized,with dark green caps,25/bag,500/case |
| | CF1500 | 15ml, Conical, Flip top cap, Sterilized (25/pkg,500/cs) |
| | CF5000 | 50ml, Conical ,Flip top cap, Sterilized (25/pkg,500/cs) |
| | DTC111500 | 50ml, Green Cap, Self-standing, Sterilized, 25/pkg, 500/case |
| | DTC1500 | 15ml, JET light green cap, Conical, Sterilized(25/pkg, 500/cs) |
| | DTC1500R | 15ml, Green Cap, Conical, Sterilized, 25/rack, 500/case |
| | DTC1550CNS | 15ml, Jet Light Green Cap, Conical, Non-Sterilized, 50/pkg, 500/case |

| | | |
|------------------------|------------|---|
| | DTC5000 | 50ml, Green Cap, Conical, Sterilized, 25/pkg, 500/case |
| | DTC5000R | 50ml, Jet Light Green Cap, Conical, Sterilized, 25/rack, 500/case |
| | DTC5050CNS | Centrifuge Tubes, 50ml, Green Cap, Conical Bottom, Non-Sterile,50/Bag,500/Case. |
| Filter Upper Cup | F0021 | 1000ml, PES, 0.22um, Sterilized (1/pkg, 24/cs). |
| | F0022 | 250ml, PES, 0.22um, Sterilized (1/pkg, 24/cs) |
| | F0025 | 500ml, PES, 0.22um, Sterilized (1/pkg, 24/cs) |
| | F0041 | 1000ml, PES, 0.45um, Sterilized (1/pkg, 24/cs) |
| | F0042 | 250ml, PES, 0.45um, Sterilized (1/pkg, 24/cs) |
| | F0045 | 500ml, PES, 0.45um, Sterilized (1/pkg, 24/cs) |
| L-Shape Cell Spreader | F4050 | Blue, Bulk Package, Sterile, Non-pyrogen (10/pkg, 500/box, 8000/cs) |
| | F4055 | Blue, Sterile, Non-pyrogen (Individually wrapped, 100/box, 1600/case) |
| Syringe Driven Filters | FE12S | 13mm, 0.22um, PES,Green, Sterilized (100/box, 800/cs) |
| | FE14S | 13mm, 0.45um, PES, Green,Sterilized (100/box,800/cs) |
| | FE32S | 30mm, 0.22um, PES,Green, Sterilized (45/box, 360/cs) |
| | FE34S | 30mm, 0.45um, PES, Green,Sterilized (45/box, 360/cs) |
| | FM32S | 30mm, 0.22um, MCE, Yellow,Sterilized (45/box, 360/cs) |
| | FM34S | 30mm, 0.45um, MCE, Yellow,Sterilized (45/box, 360/cs) |
| | FMC3022S | 30mm, 0.22um, MCE, Yellow,Sterilized (45/box, 360/cs) |
| | FMC3045S | 30mm, 0.45um, MCE, Yellow,Sterilized (45/box, 360/cs) |
| | FN14 | 13mm,0.45um,NYLON, Non-sterilized (100/bag,1000/cs) |
| | FN32S | 30mm,0.22um,NYLON,Pink, Sterilized (45/box,360/cs) |
| | FN34 | 30mm,0.45um,NYLON,Non-sterilized (50/bag,500/cs) |
| | FN34S | 30mm, 0.45um, NYLON, Pink,Sterilized (45/box, 360/cs) |
| | FNY1345 | 13mm,0.45um,NYLON, Non-sterilized (100/bag,1000/cs) |
| | FNY3022S | 30mm,0.22um,NYLON,Pink, Sterilized (45/box,360/cs) |
| | FNY3045 | 30mm,0.45um,NYLON,Non-sterilized (50/bag,500/cs) |
| | FNY3045S | 30mm, 0.45um, NYLON, Pink,Sterilized (45/box, 360/cs) |
| | FPE1322S | 13mm, 0.22um, PES,Green, Sterilized (100/box, 800/cs) |
| | FPE1345S | 13mm, 0.45um, PES, Green,Sterilized (100/box,800/cs) |
| | FPE3022S | 30mm, 0.22um, PES,Green, Sterilized (45/box, 360/cs) |
| | FPE3045S | 30mm, 0.45um, PES, Green,Sterilized (45/box, 360/cs) |

| | | |
|-----------------------|------------|--|
| | FPT2522 | 30mm, 0.22um, PTFE, Non-Sterilized (50/box, 500/cs) |
| | FPT2545 | 30mm, 0.45um, PTFE, Non-Sterilized (50/box, 500/cs) |
| | FPV1322S | 13mm, 0, 22um, PVDF,Blue, Sterilized (100/box, 800/cs) |
| | FPV3022S | 30mm, 0, 22um, PVDF, Blue,Sterilized (45/box, 360/cs) |
| | FPV3045S | 30mm, 0.45um, PVDF, Blue,Sterilized (45/box, 360/cs) |
| | FS56400-11 | Pipet,Bact,Plas,Milk,1.1mL,Multi Pack, 25/Pack,1000/Case |
| | FT22 | 30mm, 0.22um,PTFE,White, Non-Sterilized (50/bag, 500/cs). |
| | FV12S | 13mm, 0, 22um, PVDF,Blue, Sterilized (100/box, 800/cs) |
| | FV32S | 30mm, 0, 22um, PVDF, Blue,Sterilized (45/box, 360/cs) |
| | FV34S | 30mm, 0.45um, PVDF, Blue,Sterilized (45/box, 360/cs) |
| Inoculating Needles | L0000 | Innoculating Needle, White, Sterilized (25/pkg, 1000/box, 10boxes/cs) |
| | L0000-1 | Innoculating Needle, White, Sterilized, Individually wrapped, 250/box, 10boxes/cs |
| | L0000-25 | Innoculating Needle, White, Sterilized (10/pkg, 1000/box, 10boxes/cs) |
| Inoculating Loops | L0001 | 1.0ul, Blue, Sterilized (25/pkg, 1000/box, 10boxes/cs) |
| | L0010 | 10.0ul, Yellow, Sterilized (25/pkg, 1000/box, 10boxes/cs) |
| | L0010-1 | 10.0ul, Yellow, Sterilized, Individually wrapped, 250/box, 10boxes/cs |
| Open End Pipet | OEPS1 | 1ml,Open End,Sterilized (Individually wrapped,500/cs) |
| | OEPS10 | 10ml,Open End ,Sterilized (Individually wrapped ,200/cs) |
| | OEPS2 | 2ml,Open End,Sterilized (Individually wrapped,500/cs) |
| | OEPS5 | 5ml,Open End,Sterilized (Individually wrapped,200/cs) |
| Epette Pipettes | EV005 | 5ml Epette Pipette, 50 Per Inner Box, 50 Per Pack, , 500 Per Case, Sterile |
| Plastic Pasteur Pipet | P0534 | 145MM, Sterilized, Individually Wrapped, 50pcs/inner bag, 200pcs/box, 4boxes/case |
| | P0534 | 145MM, Sterilized, Individually Wrapped, 50/inner bag, 200/cs |
| | P0900 | 230MM, Sterilized, Individually Wrapped, 50pcs/inner bag, 200pcs/box, 12boxes/case |
| | P0900 | 230MM, Sterilized, Individually Wrapped, 50/inner bag, 200/cs |
| | BPS10 | Biosafe Pipette,Sero,Plas,TD,10mL,10x 1/10mL,Case of 200 |
| | BPS25 | Pipette,Sero,Plas,TD,25mL,25x 2/10mL,Case of 150 |
| | BPS5 | Biosafe Pipette, Sero, Plas, TD, 5mL, 5 x 1/10mL, Case of 200 |
| | PB01 | Pipette,Sero,Plas,TD,Multi.1mL,1x 1/100mL,Shelfpacks of 25,Case of 1000 |
| | PB1 | Pipette,Sero,Plas,TD,Multi.1mL,1x 1/10mL,Shelfpacks of 25,Case of 1000 |
| | PB10 | Pipette,Sero,Plas,TD,Multi.10mL,10x 1/10mL,Shelfpacks of 25,Case of 400 |

| | | |
|---------------------------|---|--|
| Serological Pipettes | PB10W | Pipette, Sero, Plas, TD, Multi, 10mL, 10x1/10mL, Shelfpacks of 25, Case of 400 *Wide Tip* |
| | PB2 | Pipette,Sero,Plas,TD,Multi.2mL,2x 1/100mL,Case of 1000 |
| | PB22 | Pipette,Sero,Plas,TD,25/pack,1000/cs |
| | PB25 | Pipette,Sero,Plas,TD,Multi.25mL,25x 2/10mL,Shelfpacks of 10,Case of 150 |
| | PB5 | Pipette,Sero,Plas,TD,Multi.5mL,5x 1/10mL,Shelfpacks of 25,Case of 500 |
| | PB50 | Pipette,Sero,Plas,TD,Multi.50mL,50x 5/10mL,Shelfpacks of 10,Case of 100 |
| | PM10 | Pipette, Sero, 10ml, Plas, Short, 150/Case |
| | PM25 | Pipette, Sero, Plas, 25ml Short, 100/Case |
| | PM5 | Pipette, Sero, Plas, 5ml Short, 200/Case |
| | PMB11 | Pipette,Bact,Plas,Milk,1.1mL,Shelfpacks of 25,Case of 1000 |
| | PMS11 | Pipette,Bact,Plas,Milk,1.1mL,Case of 500 |
| | PS01 | Pipette, Sero, Plas, TD, 1mL, 1x 1/100mL, Case of 500 |
| | PS1 | Pipette,Sero,Plas,TD,1mL,1x 1/10mL,Case of 500 |
| | PS10 | Pipette,Sero,Plas,TD,10mL,10x 1/10mL,Case of 200 |
| | PS100 | 100ml,Sterilized,Pink,Ind Wrap(Paper/Plastic),50/cs |
| | PS10W | 10ml,1/10ml,Sterilized,Wide tip,Orange,Ind Wrap(Paper/Plastic),200/cs |
| | PS2 | Pipette,Sero,Plas,TD,2mL,2x 1/100mL,Case of 500 |
| | PS22 | Pipette ,Sero,Plas,TD,Individually Wrapped ,500/CS |
| | PS25 | Pipette,Sero,Plas,TD,25mL,25x 2/10mL,Case of 150 |
| | PS5 | Pipette, Sero, Plas, TD, 5mL, 5 x 1/10mL, Case of 200 |
| PS50 | Pipette, Sero, Plas, TD, 50mL, 50 x 5/10mL, Case of 100 | |
| Sealing pad | S0000 | Deep well plate sealing pad, non-sterilized, 50/bag,100/cs |
| Barbed Cap | MV001 | 45mm, Non-sterilized, 1pc/bag, 200pcs/pack, 200/cs |
| Backflip Centrifuge Tubes | T0015 | Argos Backflip 1.5ml Centrifuge Tube, Non-sterilized, Dnase/Rnase free, 500/bag,5000/cs |
| | T0015S | Argos Backflip 1.5ml Centrifuge Tube, Sterilized, Dnase/Rnase free, 500/bag,5000/cs |
| 5ml Micro Centrifuge Tube | T2075A | 5mL CapSafe Microtube, Natural Color, Non-sterile, RNase, DNase & Pyrogen Free, 200/PK, 4000/case |
| | T2076A | 5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Natural, 200/PK, 4000/case |
| | T2076S-CA | 5mL Microtube, Sterilized by Gamma Radiation, RNase、DNase & Pyrogen Free, Natural, 20/PK, 200/case |
| | T2077A | 5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Blue, 200/PK, 4000/case |
| | T2078 | 5mL Non-sterile Microtube, Yellow, 250/PK, 2500/case |
| | T2078A | 5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Yellow, 200/PK, 4000/case |

| | | |
|--------------------------|-----------|---|
| | T2079A | 5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Green, 200/PK, 4000/case |
| | T2080A | 5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Rose, 200/PK, 4000/case |
| | T2081 | 5mL Non-sterile Microtube, Black, 250/PK, 2500/case |
| | T2081A | 5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Black, 200/PK, 4000/case |
| Centrifuge Tubes | SS-PH15 | 15ml, Conical, Green,Cap,Sterilized, Dnase/Rnase-free (25/pkg, 500/cs) |
| | SS-PH15R | 15ml, Conical, Green Cap,Sterilized, R/Dnase-free (25/rack, 500/cs) |
| | SS-PH50 | 50ml, Conical, Green Cap,Sterilized, Dnase/Rnase-free (25/pkg, 500/cs) |
| | SS-PH50R | 50ml, Conical, Green Cap,Sterilized, Dnase/Rnase-free (25/rack, 500/cs) |
| | TB1500 | Black 15ML Conical Tube, Sterile (50/pk, 500/cs) |
| | TB5000 | Black 50ML Conical Tube, Sterile (50/pk, 500/cs) |
| | TC1500 | 15ml,Red cap,Conical,Sterilized (25/pkg,500/cs) |
| | TC1500R | 15ml,Red cap,Conical,Sterilized (25/rack,500/cs) |
| | TC5000 | 50ml,Red cap,Conical,Sterilized (25/pkg,500/cs) |
| | TC5000R | 50ml,Red cap,Conical,Sterilized (25/rack,500/cs) |
| Pestles | V7338-901 | 0.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs |
| | V7339-901 | 1.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/case |
| Micro Centrifuge Tubes | V7811-904 | 1.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs |
| | V8185-904 | 0.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs |
| Pestle & Microtube Combo | V9950-901 | 1.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs |
| | V9951-901 | 0.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs |
| Inoculating Needles | VL0000 | Innoculating Needle, White, Sterilized (25/pkg, 1000/box, 10boxes/cs) |
| | VL0000 | Innoculating Needle, White, Sterilized (25/pkg, 1000/box, 12boxes/cs) |
| | VL0000-1 | Innoculating Needle, White, Sterilized, Individually wrapped, 250/box, 12boxes/cs |
| | VL0000-25 | Innoculating Needle, White, Sterilized (10/pkg, 1000/box, 12boxes/cs) |
| Inoculating Loops | VL0001 | 1.0ul, Blue, Sterilized (25/pkg, 1000/box, 12boxes/cs) |
| | VL0001 | 1.0ul, Blue, Sterilized (25/pkg, 1000/box, 10boxes/cs) |
| | VL0001-1 | 1.0ul, Blue, Sterilized, Individually wrapped, 250/box, 12boxes/cs |
| | VL0001-25 | 1.0ul, Blue, Sterilized (10/pkg, 1000/box, 10boxes/cs) |
| | VL0001-25 | 1.0ul, Blue, Sterilized (10/pkg, 1000/box, 12boxes/cs) |
| | VL0010 | 10.0ul, Yellow, Sterilized (25/pkg, 1000/box, 10boxes/cs) |
| | VL0010 | 10.0ul, Yellow, Sterilized (25/pkg, 1000/box, 12boxes/cs) |

| | | |
|----------------|-----------|---|
| | VL0010-1 | 10.0ul, Yellow, Sterilized, Individually wrapped, 250/box, 12boxes/cs |
| | VL0010-25 | 10.0ul, Yellow, Sterilized (10/pkg, 1000/box, 10boxes/cs) |
| | VL0010-25 | 10.0ul, Yellow, Sterilized (10/pkg, 1000/box, 12boxes/cs) |
| Roller Bottles | W2000 | Cut, No cap, No Plasma, Non-sterilized (1/pkg, 12/cs) |