



## Declaration of conformity



According to the in Vitro Diagnostic Medicine Device Directive 98/79/EC

**Manufacturer:** GUANGZHOU JET BIO-FILTRATION CO., LTD  
No.1 Doutang Road, Yonghe Economical & Technological, Development  
District, Luogang, 511356 Guangzhou, Guangdong, P.R. China  
No.173-175 Youyi Road Economical & Technological Development  
District 510730 Guangzhou, P.R. China

**Authorised**

**Representative:** TÜV Süddeutschland Group  
Ridler str.65, Munchen, 80339 Deutschland

**Description/Category:** In Vitro Diagnostic Medical Devices

**IVDD Classification :** Other Type (General)

**Products Name :** See attached list

The undersigned hereby declares that the in Vitro Diagnostic Medical Devices as specified in the attached list conform with the essential requirements listed in the Annex I of European In Vitro Diagnostic Medicine Device Directive 98/79/EC(IVDD).

This declaration of conformity is based on the European In Vitro Diagnostic Medical Device Directive 98/79/EC(IVDD). Annex III is supported by a EC-REP TUV SUD Product Service GmbH Annex V Certification. With reference to articles 1 and 3 of the IVDD.

Quality Manager

*Guangzhou*  
*Aug. 17, 2015*  
(Place and date of issue)

*York Zhang*  
(Name and signature of authorized person)



## Declaration of Conformity Directive 98/79/EC

### PRODUCTS LIST

Product Type	Cat.No.	Description and Packing
Syringe Driven Filters	23966	13mm, 0.22um, PES, Green, Non-sterile, 100/pk(box), 1000/cs
	23967	13mm, 0.45um, PES, Green, Non-sterile, 100/pk(box), 1000/cs
	23968	25mm, 0.22um, PES, Green, Non-sterile, 100/pk(box), 1000/cs
	23969	25mm, 0.45um, PES, Green, Non-sterile, 100/pk(box), 1000/cs
	23970	4mm, 0.22um, Nylon, Clear, Non-sterile, 100/pk(box), 1000/cs
	23971	4mm, 0.45um, Nylon, Clear, Non-sterile, 100/pk(box), 1000/cs
	23972	4mm, 0.22um, CA, Clear, Non-sterile, 100/pk(box), 1000/cs
	23973	4mm, 0.22um, CA, Clear, Non-sterile, 100/pk(box), 1000/cs
	23974	4mm, 0.22um, PTFE, Clear, Non-sterile, 100/pk(box), 1000/cs
	23975	4mm, 0.45um, PTFE, Clear, Non-sterile, 100/pk(box), 1000/cs
	23976	4mm, 0.22um, PVDF, Clear, Non-sterile, 100/pk(box), 1000/cs
	23977	4mm, 0.45um, PVDF, Clear, Non-sterile, 100/pk(box), 1000/cs
	23978	4mm, 0.22um, PES, Clear, Non-sterile, 100/pk(box), 1000/cs
	23979	4mm, 0.45um, PES, Clear, Non-sterile, 100/pk(box), 1000/cs
	23980	30mm, 0.22um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs
	23981	30mm, 0.45um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs
	23982	30mm, 0.22um, Cellulose, Yellow, Non-sterile, 100/pk(box), 1000/cs
	23983	30mm, 0.45um, CA, Red, Non-sterile, 100/pk(box), 1000/cs
	23984	30mm, 0.22um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs
	23985	30mm, 0.45um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs
	23986	30mm, 0.22um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs
	23987	30mm, 0.45um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs
	23988	30mm, 0.22um, PES, Green, Non-sterile, 100/pk(box), 1000/cs
	23989	30mm, 0.45um, PES, Green, Non-sterile, 100/pk(box), 1000/cs
	26142	13mm, 0.22um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs
	26143	13mm, 0.45um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs
	26144	25mm, 0.22um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs

	26145	25mm, 0.45um, PTFE, White, Non-sterile, 100/pk(box), 1000/cs
	26146	13mm, 0.22um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs
	26147	13mm, 0.45um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs
	26148	25mm, 0.22um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs
	26149	25mm, 0.45um, Nylon, Pink, Non-sterile, 100/pk(box), 1000/cs
	26150	13mm, 0.22um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs
	26151	13mm, 0.45um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs
	26152	25mm, 0.22um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs
	26153	25mm, 0.45um, PVDF, Blue, Non-sterile, 100/pk(box), 1000/cs
	26155	13mm, 0.45um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs
	26156	13mm, 0.22um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs
	26157	25mm, 0.45um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs
	26158	25mm, 0.22um, Cellulose Acetate, Red, Non-sterile, 100/pk(box), 1000/cs
Serological Pipettes	56400-11	Pipet,Bact,Plas,Milk,1.1mL,Multi Pack,25/Pack,1000/Case
	56400-22	Pipette, Sero, Plas,TD, Multi. 2.2ml, 2x 1/100mL, 25/pack, 1000/cs
	56404-11	Pipet,Bact,Plas,Milk,1.1mL,Individually Wrapped, 200/Pack,1000/Case
	56404-22	Pipette, Sero, Plas,TD, 2.2ml, 2x 1/100mL, Individually Wrapped, 500/CS
	56700-10110	Pipet,Sero,Plas,TD,Canister Pack,10ml,Plugged,Sterile,50/Canister,500/Case
	56700-11100	Pipet,Sero,Plas,TD,Canister Pack,1ml,Plugged,Sterile,50/Canister,500/Case
	56700-21100	Pipet,Sero,Plas,TD,Canister,2ml,Plugged,50/Canister,500/Case
	56700-5110	Pipet,Sero,Plas,TD,Canister,5ml,Plugged,50/Canister,500/Case
	56800-10110	Pipette,Sero,Plas,TD,Multi.10mL,10x 1/10mL,Shelfpacks of 25,Case of 500
	56800-1110	Pipette,Sero,Plas,TD,Multi.1mL,1x 1/10mL,Shelfpacks of 25,Case of 1000
	56800-11100	Pipet, Sero, Plas, TD, Multi, 1ml, Plugged, Sterile, 25/Pack, 1000/Case
	56800-21100	Pipet,Sero,Plas,TD,Multi Pack,2ml, 25/Pack,500/Case
	56800-25210	Pipet,Sero,Plas,TD,Multi Pack,25mL,Plugged, Sterile,10/pack,200/case.
	56800-50510	Pipette,Sero,Plas,TD,Multi.50mL,50x 5/10mL,Shelfpacks of 10,Case of 100
	56800-5110	Pipet,Sero,Plas,TD,Multi Pack,5ml,Plastic, Plugged, Sterile, 25/Pack,500/Case
	56825-10110	Pipet, Sero, Plas, TD, Multi Pack, 10mL,Plugged, Sterile, 25/Pack, 500/Case *Wide Tip*
	56900-10010	100ml,Sterilized,Pink,Ind Wrap(Paper/Plastic),50/cs
	56900-10110	Pipette,Sero,Plas,TD,10ml,Individually Wrapped, Plugged, Sterile,200/Case (Shelfpacks of 50)

	56900-1110	Pipet,Sero,Plas,TD,1mL,Individually Wrapped, Plugged, Sterile,1000/Case (Shelfpacks of200)
	56900-11100	Pipet,Sero,Plas,TD,1ml,Individually Wrapped, Plugged, Sterile,1000/Case (Shelfpacks of 200)
	56900-21100	Pipet,Sero,Plas,TD,2ml,Individually Wrapped, Plugged, Sterile,500/Case (Shelfpacks of 125)
	56900-25210	Pipette,Sero,Plas,TD,25mL,Individually Wrapped, Plugged, Sterile,200/Case (Shelfpacks of 50)
	56900-50510	Pipet, Sero, Plas, TD, 50ml,Individually Wrapped, Plugged,Sterile, 100/Case (Shelfpacks of 25)
	56900-5110	Pipet,Sero,Plas,TD,5mL,Individually Wrapped, Plugged, Sterile,200/Case (Shelfpacks of 50)
Aspirating Pipets	APS10	10ml,Sterilized ( Individually wrapped ,200/cs)
	APS2	2ml,Sterilized (Individually wrapped,500/cs)
	APS5	5ml,Sterilized (Individually wrapped,200/cs)
Reservoir Bottle	B0001	1000ml, Non-sterilized (1/pkg, 24/cs)
Pipetting Reservoir	B3110-100	100mL, Sterile,DNA-RNase & pyrogen free, 5/pk, 100pcs/case.
	B3110-100NS	100mL, Bulk Pack , Non Sterile, DNA-RNase & pyrogen free, 100pcs/bag, 1bag/case.
	B3110-50	100mL, Sterile, Individually Wrapped, DNA-RNase & pyrogen free, 50pcs/case.
	B3125-100	25mL,Sterile,DNA-RNase & pyrogen free,5/pk, 100/case.
	B3125-100NS	25mL , Non Sterile, DNA-RNase & pyrogen free, Bulk Pack, 100pcs/bag, 1bag/case.
	B3125-50	25mL, Sterile, DNA-RNase & pyrogen free, Individually Wrapped, 50pcs/case.
	B3150-100	50mL, Sterile,DNA-RNase & pyrogen free, 5/pk, 100pcs/case.
	B3150-100NS	50mL, Bulk Pack , Non Sterile, DNA-RNase & pyrogen free, 100pcs/bag, 1bag/case.
	B3150-50	50mL, Sterile, DNA-RNase & pyrogen free, Individually Wrapped, 50pcs/case.
12 Channel basin	B3134	12 Channel non-sterile basin, individually packed, 50/cs
	B3135	12 Channel sterile basin, Individually packed, 20/box, 12boxes/case
Vacuum Filtration	BMC2225	250ml, MCE, 0.22um, Sterilized (1/pkg, 12/cs)
	BMC2250	500ml, MCE, 0.22um, Sterilized (1/pkg, 12/cs)
	BMC4525	250ml, MCE, 0.45um, Sterilized (1/pkg, 12/cs)
	BMC4550	500ml, MCE, 0.45um, Sterilized (1/pkg, 12/cs)
	BNY2225	250ml,NYLON,0.22um, Sterilized (1/pkg,12/cs)
	BNY2250	500ml,NYLON,0.22um, Sterilized (1/pkg,12/cs)
	BNY4525	250ml, NYLON, 0.45um, Sterilized (1/pkg, 12/cs)
	BNY4550	500ml, NYLON, 0.45um, Sterilized (1/pkg, 12/cs)
	BPE2210	1000ml, PES, 0.22um, Sterilized (1/pkg, 12/cs)
	BPE2215	150ml, PES, 0.22um, Sterilized (1/pkg, 12/cs).

	BPE2225	250ml, PES, 0.22um, Sterilized (1/pkg,12/cs)
	BPE2250	500ml, PES, 0.22um, Sterilized (1/pkg, 12/cs)
	BPE4510	1000ml, PES, 0.45um, Sterilized (1/pkg, 12/cs)
	BPE4515	150ml, PES, 0.45um, Sterilized (1/pkg, 12/cs)
	BPE4525	250ml, PES, 0.45um, Sterilized (1/pkg, 12/cs)
	BPE4550	500ml, PES, 0.45um, Sterilized (1/pkg, 12/cs)
Vacuum Filtration	BPV2210	1000ml, PVDF, 0.22um, Sterilized (1/pkg, 12/cs)
	BPV2225	250ml, PVDF, 0.22um, Sterilized (1/pkg, 12/cs)
	BPV2250	500ml, PVDF, 0.22um, Sterilized (1/pkg, 12/cs)
	BPV4510	1000ml, PVDF, 0.45um, Sterilized (1/pkg, 12/cs)
	BPV4525	250ml, PVDF, 0.45um, Sterilized (1/pkg, 12/cs)
	BPV4550	500ml, PVDF, 0.45um, Sterilized (1/pkg, 12/cs)
Cell Scraper	C2500	25cm ,Cell Scraper,Purple,Sterilized,Non-pyrogenic,Lifter Blade( Individually wrapped, 100pcs per bag/case)
	C3900	39cm ,Cell Scraper,Sterilized,Non-pyrogenic,Lifter Blade( Individually wrapped, 100pcs per
Cell Strainer	CS040	40µm Cell Strainer, Blue, Sterile, Individually Packaged, 50/Case
	CS070	70µm Cell Strainer, White, Sterile, Individually Packaged, 50/Case
	CS100	100µm Cell Strainer, Yellow, Sterile, Individually Packaged, 50/Case
Pipetting Reservoir	D3110-100	Denville Basins, 100mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free (100/case). Denville # 1005110.
	D3110-100NS	Denville Basins, 100mL, Bulk Pack, Non Sterile, 100/cs. Denville # 1005118
	D3110-50	Denville Basins, 100mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/Case. Denville#1005108.
	D3125-100	Denville Basins, 25mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free 100/case. Denville # 1005120
	D3125-100NS	Denville Basins, 25mL, Bulk Pack, Non-Sterile, 100/cs. Denville # 1005121
	D3125-50	Denville Basins, 25mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/case. Denville # 1005119
	D3150-100	Denville Basins, 50mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free 100/case. Denville # 1005109
	D3150-100NS	Denville Basins, 50mL, Bulk Pack, Non Sterile, 100/cs. Denville # 1005111
	D3150-50	Denville Basins, 50mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/case. Denville # 1005122
	DB3110-50	ClearLine Basins, 100mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free 50/Case
	DB3125-100	ClearLine Basins, 25mL, 5/pk, Sterile, DNA-RNase & Pyrogen Free 100/case
	DB3125-100NS	ClearLine Basins, 25mL, Bulk Pack, Non-Sterile, 100/cs
	DB3125-50	ClearLine Basins, 25mL, Sterile, Individually, Wrapped, DNA-RNase & Pyrogen Free 50/case

	DB3150-50	ClearLine Basins, 50mL, Sterile, Individually Wrapped, DNA-RNase & Pyrogen Free, 50/case
Cell Scraper	DC2500	25cm ,Cell Scraper, Purple, Sterilized, Nonpyrogenic, Lifter Blade (Individually wrapped, 100pcs per bag/case)
	DC3900	39cm ,Cell Scraper, Purple, Sterilized, Nonpyrogenic, Lifter Blade (Individually wrapped, 100pcs per case)
Tissue Culture Dish	DDT1000	Tissue culture dishes, 100mm, surface treated, sterile, 10/bag, 300/case
	DDT3500	Tissue culture dishes, 35mm, surface treated, sterile, 10/bag, 960/case
	DDT7000	Tissue culture dishes, 70mm, surface treated, sterile, 10/bag, 600/case
Tissue Culture Flask	DF1105	Tissue culture flask, 50ml, 25cm <sup>2</sup> , non-vented, sterile with treated surface,10/bag 200/case
	DF1125	Tissue culture flask, 250ml, 75cm <sup>2</sup> ,non-vented, sterile with treated surface,5/bag ,100/case
	DF1160	Tissue culture flask, 600ml, 182cm <sup>2</sup> ,non-vented, sterile with treated surface,5/bag, 40/case
	DF1205	Tissue Culture Flask, 50ml with Vented Cap, Sterile with Treated Surface,10/Bag, 200/Case
	DF1225	Tissue Culture Flask, 250ml, with Vented Cap, Sterilize with Treated Surface, 5/Bag,100/Case.
	DF1260	Tissue Culture Flask, 600ml with Vented Cap, Sterile with Treated Surface, 5/Bag 40/Case
	F1205	50ml, Surface-treated, w/vent cap, Sterilized (10/pkg, 200/cs)
Tissue Culture Plate	DP012096	96 well, U-bottom, Surface-treated, Sterilized, 1/pkg, 100/cases
	DP1006	6well, Surface-treated, Sterilized (1/pkg, 100/cs)
	DP1012	12 well, Surface-treated, Sterilized, 1/pkg, 100/cs
	DP1024	24 well, Surface-treated, Sterilized, 1/pkg, 100/case
	DP1096	96 well, Surface-treated, Sterilized, 1/pkg, 100/case
	P0096	96well, Non-treated, Sterilized (1/pkg, 100/cs)
	P1006	6well, Surface-treated, Sterilized (1/pkg, 100/cs)
	P1012	12well, Surface-treated, Sterilized (1/pkg, 100/cs)
	P1024	24well, Surface-treated, Sterilized (1/pkg, 100/cs)
	P1096	96well, Surface-treated, Sterilized (1/pkg, 100/cs)
	P2096	96well,without lid, non-treated and non-steriled,(10/bag,100/cs).
Centrifuge Tubes	4503	50ml,Conical,Sterilized,with dark green caps,25/bag,500/case
	CF1500	15ml, Conical, Flip top cap, Sterilized (25/pkg,500/cs)
	CF5000	50ml, Conical ,Flip top cap, Sterilized (25/pkg,500/cs)
	DTC111500	50ml, Green Cap, Self-standing, Sterilized, 25/pkg, 500/case
	DTC1500	15ml, JET light green cap, Conical, Sterilized(25/pkg, 500/cs)
	DTC1500R	15ml, Green Cap, Conical, Sterilized, 25/rack, 500/case
	DTC1550CNS	15ml, Jet Light Green Cap, Conical, Non-Sterilized, 50/pkg, 500/case

	DTC5000	50ml, Green Cap, Conical, Sterilized, 25/pkg, 500/case
	DTC5000R	50ml, Jet Light Green Cap, Conical, Sterilized, 25/rack, 500/case
	DTC5050CNS	Centrifuge Tubes, 50ml, Green Cap, Conical Bottom, Non-Sterile,50/Bag,500/Case.
Filter Upper Cup	F0021	1000ml, PES, 0.22um, Sterilized (1/pkg, 24/cs).
	F0022	250ml, PES, 0.22um, Sterilized (1/pkg, 24/cs)
	F0025	500ml, PES, 0.22um, Sterilized (1/pkg, 24/cs)
	F0041	1000ml, PES, 0.45um, Sterilized (1/pkg, 24/cs)
	F0042	250ml, PES, 0.45um, Sterilized (1/pkg, 24/cs)
	F0045	500ml, PES, 0.45um, Sterilized (1/pkg, 24/cs)
L-Shape Cell Spreader	F4050	Blue, Bulk Package, Sterile, Non-pyrongen (10/pkg, 500/box, 8000/cs)
	F4055	Blue, Sterile, Non-pyrongen (Individually wrapped, 100/box, 1600/case)
Syringe Driven Filters	FE12S	13mm, 0.22um, PES,Green, Sterilized (100/box, 800/cs)
	FE14S	13mm, 0.45um, PES, Green,Sterilized (100/box,800/cs)
	FE32S	30mm, 0.22um, PES,Green, Sterilized (45/box, 360/cs)
	FE34S	30mm, 0.45um, PES, Green,Sterilized (45/box, 360/cs)
	FM32S	30mm, 0.22um, MCE, Yellow,Sterilized (45/box, 360/cs)
	FM34S	30mm, 0.45um, MCE, Yellow,Sterilized (45/box, 360/cs)
	FMC3022S	30mm, 0.22um, MCE, Yellow,Sterilized (45/box, 360/cs)
	FMC3045S	30mm, 0.45um, MCE, Yellow,Sterilized (45/box, 360/cs)
	FN14	13mm,0.45um,NYLON, Non-sterilized (100/bag,1000/cs)
	FN32S	30mm,0.22um,NYLON,Pink, Sterilized (45/box,360/cs)
	FN34	30mm,0.45um,NYLON,Non-sterilized (50/bag,500/cs)
	FN34S	30mm, 0.45um, NYLON, Pink,Sterilized (45/box, 360/cs)
	FNY1345	13mm,0.45um,NYLON, Non-sterilized (100/bag,1000/cs)
	FNY3022S	30mm,0.22um,NYLON,Pink, Sterilized (45/box,360/cs)
	FNY3045	30mm,0.45um,NYLON,Non-sterilized (50/bag,500/cs)
	FNY3045S	30mm, 0.45um, NYLON, Pink,Sterilized (45/box, 360/cs)
	FPE1322S	13mm, 0.22um, PES,Green, Sterilized (100/box, 800/cs)
	FPE1345S	13mm, 0.45um, PES, Green,Sterilized (100/box,800/cs)
	FPE3022S	30mm, 0.22um, PES,Green, Sterilized (45/box, 360/cs)
	FPE3045S	30mm, 0.45um, PES, Green,Sterilized (45/box, 360/cs)

	FPT2522	30mm, 0.22um, PTFE, Non-Sterilized (50/box, 500/cs)
	FPT2545	30mm, 0.45um, PTFE, Non-Sterilized (50/box, 500/cs)
	FPV1322S	13mm, 0, 22um, PVDF,Blue, Sterilized (100/box, 800/cs)
	FPV3022S	30mm, 0, 22um, PVDF, Blue,Sterilized (45/box, 360/cs)
	FPV3045S	30mm, 0.45um, PVDF, Blue,Sterilized (45/box, 360/cs)
	FS56400-11	Pipet,Bact,Plas,Milk,1.1mL,Multi Pack, 25/Pack,1000/Case
	FT22	30mm, 0.22um,PTFE,White, Non-Sterilized (50/bag, 500/cs).
	FV12S	13mm, 0, 22um, PVDF,Blue, Sterilized (100/box, 800/cs)
	FV32S	30mm, 0, 22um, PVDF, Blue,Sterilized (45/box, 360/cs)
	FV34S	30mm, 0.45um, PVDF, Blue,Sterilized (45/box, 360/cs)
Inoculating Needles	L0000	Innoculating Needle, White, Sterilized (25/pkg, 1000/box, 10boxes/cs)
	L0000-1	Innoculating Needle, White, Sterilized, Individually wrapped, 250/box, 10boxes/cs
	L0000-25	Innoculating Needle, White, Sterilized (10/pkg, 1000/box, 10boxes/cs)
Inoculating Loops	L0001	1.0ul, Blue, Sterilized (25/pkg, 1000/box, 10boxes/cs)
	L0010	10.0ul, Yellow, Sterilized (25/pkg, 1000/box, 10boxes/cs)
	L0010-1	10.0ul, Yellow, Sterilized, Individually wrapped, 250/box, 10boxes/cs
Open End Pipet	OEPS1	1ml,Open End,Sterilized (Individually wrapped,500/cs)
	OEPS10	10ml,Open End ,Sterilized (Individually wrapped ,200/cs)
	OEPS2	2ml,Open End,Sterilized (Individually wrapped,500/cs)
	OEPS5	5ml,Open End,Sterilized (Individually wrapped,200/cs)
Epette Pipettes	EV005	5ml Epette Pipette, 50 Per Inner Box, 50 Per Pack, , 500 Per Case, Sterile
Plastic Pasteur Pipet	P0534	145MM, Sterilized, Individually Wrapped, 50pcs/inner bag, 200pcs/box, 4boxes/case
	P0534	145MM, Sterilized, Individually Wrapped, 50/inner bag, 200/cs
	P0900	230MM, Sterilized, Individually Wrapped, 50pcs/inner bag, 200pcs/box, 12boxes/case
	P0900	230MM, Sterilized, Individually Wrapped, 50/inner bag, 200/cs
	BPS10	Biosafe Pipette,Sero,Plas,TD,10mL,10x 1/10mL,Case of 200
	BPS25	Pipette,Sero,Plas,TD,25mL,25x 2/10mL,Case of 150
	BPS5	Biosafe Pipette, Sero, Plas, TD, 5mL, 5 x 1/10mL, Case of 200
	PB01	Pipette,Sero,Plas,TD,Multi.1mL,1x 1/100mL,Shelfpacks of 25,Case of 1000
	PB1	Pipette,Sero,Plas,TD,Multi.1mL,1x 1/10mL,Shelfpacks of 25,Case of 1000
	PB10	Pipette,Sero,Plas,TD,Multi.10mL,10x 1/10mL,Shelfpacks of 25,Case of 400



Serological Pipettes	PB10W	Pipette, Sero, Plas, TD, Multi, 10mL, 10x1/10mL, Shelfpacks of 25, Case of 400 *Wide Tip*
	PB2	Pipette,Sero,Plas,TD,Multi.2mL,2x 1/100mL,Case of 1000
	PB22	Pipette,Sero,Plas,TD,25/pack,1000/cs
	PB25	Pipette,Sero,Plas,TD,Multi.25mL,25x 2/10mL,Shelfpacks of 10,Case of 150
	PB5	Pipette,Sero,Plas,TD,Multi.5mL,5x 1/10mL,Shelfpacks of 25,Case of 500
	PB50	Pipette,Sero,Plas,TD,Multi.50mL,50x 5/10mL,Shelfpacks of 10,Case of 100
	PM10	Pipette, Sero, 10ml, Plas, Short, 150/Case
	PM25	Pipette, Sero, Plas, 25ml Short, 100/Case
	PM5	Pipette, Sero, Plas, 5ml Short, 200/Case
	PMB11	Pipette,Bact,Plas,Milk,1.1mL,Shelfpacks of 25,Case of 1000
	PMS11	Pipette,Bact,Plas,Milk,1.1mL,Case of 500
	PS01	Pipette, Sero, Plas, TD, 1mL, 1x 1/100mL, Case of 500
	PS1	Pipette,Sero,Plas,TD,1mL,1x 1/10mL,Case of 500
	PS10	Pipette,Sero,Plas,TD,10mL,10x 1/10mL,Case of 200
	PS100	100ml,Sterilized,Pink,Ind Wrap(Paper/Plastic),50/cs
	PS10W	10ml,1/10ml,Sterilized,Wide tip,Orange,Ind Wrap(Paper/Plastic),200/cs
	PS2	Pipette,Sero,Plas,TD,2mL,2x 1/100mL,Case of 500
	PS22	Pipette ,Sero,Plas,TD,Individually Wrapped ,500/CS
	PS25	Pipette,Sero,Plas,TD,25mL,25x 2/10mL,Case of 150
	PS5	Pipette, Sero, Plas, TD, 5mL, 5 x 1/10mL, Case of 200
PS50	Pipette, Sero, Plas, TD, 50mL, 50 x 5/10mL, Case of 100	
Sealing pad	S0000	Deep well plate sealing pad, non-sterilized, 50/bag,100/cs
Barbed Cap	MV001	45mm, Non-sterilized, 1pc/bag, 200pcs/pack, 200/cs
Backflip Centrifuge Tubes	T0015	Argos Backflip 1.5ml Centrifuge Tube, Non-sterilized, Dnase/Rnase free, 500/bag,5000/cs
	T0015S	Argos Backflip 1.5ml Centrifuge Tube, Sterilized, Dnase/Rnase free, 500/bag,5000/cs
5ml Micro Centrifuge Tube	T2075A	5mL CapSafe Microtube, Natural Color, Non-sterile, RNase, DNase & Pyrogen Free, 200/PK, 4000/case
	T2076A	5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Natural, 200/PK, 4000/case
	T2076S-CA	5mL Microtube, Sterilized by Gamma Radiation, RNase、DNase & Pyrogen Free, Natural, 20/PK, 200/case
	T2077A	5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Blue, 200/PK, 4000/case
	T2078	5mL Non-sterile Microtube, Yellow, 250/PK, 2500/case
	T2078A	5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Yellow, 200/PK, 4000/case

	T2079A	5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Green, 200/PK, 4000/case
	T2080A	5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Rose, 200/PK, 4000/case
	T2081	5mL Non-sterile Microtube, Black, 250/PK, 2500/case
	T2081A	5mL Non-sterile Microtube, RNase、DNase & Pyrogen Free, Black, 200/PK, 4000/case
Centrifuge Tubes	SS-PH15	15ml, Conical, Green,Cap,Sterilized, Dnase/Rnase-free (25/pkg, 500/cs)
	SS-PH15R	15ml, Conical, Green Cap,Sterilized, R/Dnase-free (25/rack, 500/cs)
	SS-PH50	50ml, Conical, Green Cap,Sterilized, Dnase/Rnase-free (25/pkg, 500/cs)
	SS-PH50R	50ml, Conical, Green Cap,Sterilized, Dnase/Rnase-free (25/rack, 500/cs)
	TB1500	Black 15ML Conical Tube, Sterile (50/pk, 500/cs)
	TB5000	Black 50ML Conical Tube, Sterile (50/pk, 500/cs)
	TC1500	15ml,Red cap,Conical,Sterilized (25/pkg,500/cs)
	TC1500R	15ml,Red cap,Conical,Sterilized (25/rack,500/cs)
	TC5000	50ml,Red cap,Conical,Sterilized (25/pkg,500/cs)
	TC5000R	50ml,Red cap,Conical,Sterilized (25/rack,500/cs)
Pestles	V7338-901	0.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs
	V7339-901	1.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/case
Micro Centrifuge Tubes	V7811-904	1.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs
	V8185-904	0.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs
Pestle & Microtube Combo	V9950-901	1.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs
	V9951-901	0.5ml, Sterilized, DNase/RNase free, individually packed, 100/bag, 10bags/cs
Inoculating Needles	VL0000	Innoculating Needle, White, Sterilized (25/pkg, 1000/box, 10boxes/cs)
	VL0000	Innoculating Needle, White, Sterilized (25/pkg, 1000/box, 12boxes/cs)
	VL0000-1	Innoculating Needle, White, Sterilized, Individually wrapped, 250/box, 12boxes/cs
	VL0000-25	Innoculating Needle, White, Sterilized (10/pkg, 1000/box, 12boxes/cs)
Inoculating Loops	VL0001	1.0ul, Blue, Sterilized (25/pkg, 1000/box, 12boxes/cs)
	VL0001	1.0ul, Blue, Sterilized (25/pkg, 1000/box, 10boxes/cs)
	VL0001-1	1.0ul, Blue, Sterilized, Individually wrapped, 250/box, 12boxes/cs
	VL0001-25	1.0ul, Blue, Sterilized (10/pkg, 1000/box, 10boxes/cs)
	VL0001-25	1.0ul, Blue, Sterilized (10/pkg, 1000/box, 12boxes/cs)
	VL0010	10.0ul, Yellow, Sterilized (25/pkg, 1000/box, 10boxes/cs)
	VL0010	10.0ul, Yellow, Sterilized (25/pkg, 1000/box, 12boxes/cs)

	VL0010-1	10.0ul, Yellow, Sterilized, Individually wrapped, 250/box, 12boxes/cs
	VL0010-25	10.0ul, Yellow, Sterilized (10/pkg, 1000/box, 10boxes/cs)
	VL0010-25	10.0ul, Yellow, Sterilized (10/pkg, 1000/box, 12boxes/cs)
Roller Bottles	W2000	Cut, No cap, No Plasma, Non-sterilized (1/pkg, 12/cs)