

## TECHNICAL DATA SHEET

---

<b>Part number</b>	141936
<b>Description</b>	Microtiter Drawer Racks



### Description :

- ✓ Material : stainless steel
- ✓ For Microtiter Plates up to 86x128x26 mm
- ✓ Open Design
- ✓ Grip Rail
- ✓ Without Safety Stop
- ✓ Base of drawer closed
- ✓ 7D/9H (7 deep/ 9 high)
- ✓ External dimensions (WxDxH): 138x618x276 mm

### Piece by sales :

1 Rack

---