

NEW




EPPi® 92

EPPi® 67

3 x 3 holes for **25 ml**
Eppendorf Conical Tubes*

5 x 5 holes for **5 ml**
Eppendorf Conical Tubes*

Hole diameter 30 mm as in EPPi® Boxes for 50 ml centrifuge tubes. External dimensions 133 x 133 x 92 mm.



Hole diameter 17 mm as in EPPi® Boxes for 15 ml centrifuge tubes. External dimensions 133 x 133 x 67 mm.



Cat. No.	Colour
KES92-CO9NA	natural
KES92-CO9NY	neon-yellow
KES92-CO9NB	neon-blue
KES92-CO9NG	neon-green
KES92-CO9NP	neon-red/pink
KES92-CO9NO	neon-orange

Other colours on request.

Cat. No.	Colour
KES67-CO25NA	natural
KES67-CO25NY	neon-yellow
KES67-CO25NB	neon-blue
KES67-CO25NG	neon-green
KES67-CO25NP	neon-red/pink
KES67-CO25NO	neon-orange

Other colours on request.

Made of polypropylene (PP) with two hole inserts.
With numeric ID code. Alphanumeric ID code on request.

* Eppendorf® and Eppendorf Tubes® are registered trademarks of Eppendorf AG Germany.



EPPi[®] Box made of Polypropylene for 15 ml and 50 ml centrifuge tubes



EPPi[®] 128 / 5x5 • height 128 mm fix

Inner height 125 mm. With two hole inserts 5x5 holes for 25 centrifuge tubes 15 ml (up to 17 mm Ø).
With numeric ID code. Alphanumeric ID code on request.

Cat. No.	Colour
KES128-CO25NA	natural
KES128-CO25NY	neon-yellow
KES128-CO25NB	neon-blue
KES128-CO25NG	neon-green
KES128-CO25NP	neon-red/pink
KES128-CO25NO	neon-orange
KES128-CO25WH	white

Cat. No.	Colour
KES128-CO25YE	yellow
KES128-CO25BU	blue
KES128-CO25GR	green
KES128-CO25RE	red
KES128-CO25VT	violet
KES128-CO25BL	black
KES128-CO25BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPPi[®] 128 / 3x3 • height 128 mm fix

Inner height 125 mm. With two hole inserts 3x3 holes for 9 centrifuge tubes 50 ml (up to 30 mm Ø).
With numeric ID code. Alphanumeric ID code on request.

Cat. No.	Colour
KES128-CO9NA	natural
KES128-CO9NY	neon-yellow
KES128-CO9NB	neon-blue
KES128-CO9NG	neon-green
KES128-CO9NP	neon-red/pink
KES128-CO9NO	neon-orange
KES128-CO9WH	white

Cat. No.	Colour
KES128-CO9YE	yellow
KES128-CO9BU	blue
KES128-CO9GR	green
KES128-CO9RE	red
KES128-CO9VT	violet
KES128-CO9BL	black
KES128-CO9BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPi[®] 128 / 10 L • height 128 mm fix

Inner height 125 mm. With two hole inserts for 10 centrifuge tubes 50 ml (up to 30 mm Ø).
Without ID code.

Cat. No.	Colour
KEA128-CO10NA	natural
KEA128-CO10NY	neon-yellow
KEA128-CO10NB	neon-blue
KEA128-CO10NG	neon-green
KEA128-CO10NP	neon-red/pink
KEA128-CO10NO	neon-orange
KEA128-CO10WH	white

Cat. No.	Colour
KEA128-CO10YE	yellow
KEA128-CO10BU	blue
KEA128-CO10GR	green
KEA128-CO10RE	red
KEA128-CO10VT	violet
KEA128-CO10BL	black
KEA128-CO10BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching stainless steel racks on request.

