

# Microtiter Drawer Racks

For MT-Plates 86 x 128 mm (WxD) • Stainless steel

For Microtiter Plates or DeepWells of all familiar brands.  
Base of drawers closed. Draw front closed to maximum.



## Up to 26 mm plate height • with grip rail Open Design

Cat. No. SM26 - x PC V	Plate Dimensions							Without Stop	
	3x13	4x13	5x13	6x13	7x13	8x13	O	M	
13							398	444	
12							368	410	
11							337	376	
10							307	342	
09							276	308	
08							246	274	
07							215	240	
06							185	206	
137									
	O = Without Stop	274	360	446	532	618	704		
	M = With	294	380	466	552	638	724		

Example: Cat. No. SM26 - 5 x 1 3 PCOV

MT drawer rack for 65 microtiter plates **without** safety stop.  
Dimensions: WxDxH = 137x446x398 mm



## Up to 33 mm plate height • with hand grip Open Design

Cat. No. SM33 - x PB V	Plate Dimensions							Without Stop	
	3x10	4x10	5x10	6x10	7x10	8x10	O	M	
10							377	412	
09							339	371	
08							302	330	
07							264	289	
06							227	248	
05							189	207	
04							152	166	
137									
	O = Without Stop	294	380	466	552	638	724		
	M = With	294	380	466	552	638	724		

## Up to 39 mm plate height • with hand grip Open Design

Cat. No. SM39 - x PB V	Plate Dimensions							Without Stop	
	3x09	4x09	5x09	6x09	7x09	8x09	O	M	
09							393	425	
08							350	378	
07							306	331	
06							263	284	
05							219	237	
04							176	190	
137									
	O = Without Stop	294	380	466	552	638	724		
	M = With	294	380	466	552	638	724		

