



sharpsafe
SHARPSAFE 5^{ÈME} GÉNÉRATION 

PRODUCT DATA SHEET

General data

| | |
|------------------------------|--|
| Item code | SHA-51902610 |
| GTIN code | 5056271508443 |
| Product Name | Sharpsafe® 5 th GEN. Quiver 7,5 liters |
| Description | Sharps container 7,5 liters componed by 1 body + 1 lid with hinges |
| Brand | Sharpsafe® |
| Business Unit | Safety Solutions |
| Category of safety packaging | Sharps container |
| Indication | Segregation, protection, transport, disposal of solid sharps biohazard waste |

Characteristics data

| | |
|--------------------------------|--|
| Nominal volume | 7,50 l |
| Total volume | 7,30 l |
| Usefull volume | 6,20 l |
| Maximum weight | 2 Kg |
| Body background aspect | Neutral dark grey |
| Lid color | Yellow |
| Label or printing color | Yellow |
| Inner bag color | NC |
| Recognizable dominant color | Yellow |
| Marking languages | FR/NL |
| Dimensions (assembled) in mm : | |
| Long | 150 mm |
| Wide | 100 mm |
| Height | 610 mm |
| Opening | 140 x 90 mm |
| Fill line height | 502 mm |
| Weight | 880 g |
| Material | Polypropylene recycled 45% (body) + Polypropylene virgin (55% body & 100% lid) |
| Phtalates & Latex free | Yes |
| Carbon footprint reduction | 40% |
| Autoclavable | Yes |
| Autoclavability conditions | Not autoclavable |
| Temporary closure | Yes |
| Permanent closure | Yes |
| Handles | 1 |
| Availability Accessories | Stands, trolley brackets (see brochure and compatibility tables) |



Logistic data

| | |
|---------------------------|--------------------|
| Carton quantity | 20 |
| Carton weight | 17 Kg |
| Carton dimensions (LxWxH) | 765 x 450 x 615 mm |
| Carton volume | 211,7138 dm3 |
| Palet quantity | 120 units |
| Pallet weight | 117 Kg |
| Pallet dimension (LxWxH) | 120 x 80 x 230 cm |

Certifications data

| | |
|-----------------------------|--|
| ADR (UN3291) | Yes |
| EN ISO23907-1:2019 | Yes |
| Quality mark | BSI Kitemark - (NF 302 in progress) |
| TRBA250 | Yes |
| NFX30-507:2018 | NC |
| CE MARK | 94/62/EC packaging and packaging waste |
| Quality management system : | |
| ISO 9001 | Yes |
| ISO 13485 | Yes |
| ISO 14001 | Yes |
| ISO 45001 | Yes |

Manufacture data

| | |
|-------------|-----------------------|
| Made in | United Kingdom |
| By | Frontier Plastics Ltd |
| Marketed in | A/B/CH/D/F/LUX/NL |
| By | Hospidex Group |