




sharpsafe
SHARPSAFE
5^{ÈME} GENERATION 

PRODUCT DATA SHEET

General data	
Item code	SHA-51952610 (without funnel) - SHA-51951610 (with funnel)
GTIN code	5056271508504 (without funnel) - 5056271508764 (with funnel)
Product Name	Sharpsafe® 5 th GEN. Quiver XL 25 liters (funnel in option)
Description	Sharps container 25 liters composed by 1 body + 1 lid with hinges
Brand	Sharpsafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste
Characteristics data	
Nominal volume	25 l
Total volume	25,20 l
Usefull volume	19,60 l
Maximum weight	9 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
Dimensions (assembled) in mm :	
Long	310 mm
Wide	230 mm
Height	790 mm
Opening	260 x 165 mm
Fill line height	669 mm
Weight	1 484 g (1 604 g with funnel)
Material	Polypropylene recycled 30% (body) + Polypropylene virgin (70% body & 100% lid)
Phtalates & Latex free	Yes
Carbon footprint reduction	30%
Autoclavable	Yes
Autoclavability conditions	Not autoclavable
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availability Accessories	Stands, trolley brackets (see brochure and compatibility tables)
	
Logistic data	
Carton quantity	10
Carton weight	15,10 Kg
Carton dimensions (LxWxH)	455 x 340 x 914 mm
Carton volume	141,3958 dm3
Palet quantity	100 units
Pallet weight	166 Kg
Pallet dimension (LxWxH)	120 x 80 x 230 cm
Certifications data	
ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark - (NF 302 in progress)
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes
Manufacture data	
Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group