

1. Part Number: 8049H-2000

2. Photograph:



3. Product Information:

Description:	SteriWare® Jug 2000ml
Use:	Single use, disposable jug
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Internal Diameter:	129mm
Max. Diameter of Base:	157mm
Overall Height:	167mm
Nominal Weight:	232g
Moulding & Packing Environment:	Moulded & packed in a Class 100,000 cleanroom.
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (minimum dose 25kGy)
BSE/TSE Free	Yes

Number of jugs per Box: 20
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This sampler has not been ATEX tested – Contact Sampling Systems for further information

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Sep-19