

**1. Part Number:** 8049H-2000S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Jug 2000ml - sterile
Use:	Single use, disposable jug
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Internal Diameter:	129mm
Max. Diameter of Base:	157mm
Overall Height:	167mm
Nominal Weight:	232g
Moulding & Packing Environment:	Moulded & packed in a Class 100,000 cleanroom.
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (minimum dose 25kGy)
BSE/TSE Free	Yes

Number of Samplers per Box: 16  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years from the date of manufacture  
ATEX: This sampler has not been ATEX tested – Contact Sampling Systems for further information

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Sep-19