

1. Part Numbers: A125-2000

2. Photograph:




3. Product Information:

Description	PTFE Bottle 2000ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	2000ml
Overall Diameter:	120mm
Internal Diameter of Neck:	60mm
Overall Height:	240mm
Nominal Weight:	895g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com