

1. Part Number: 8890P-500

2. Photograph:



3. Product Information:

| | |
|---|--|
| Description: | SteriWare® Powder Bottle 500ml |
| Use: | Single use, disposable bottle |
| Material of Construction: | Polypropylene (PP) – opaque <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulation 10/2011 • Conforms to EC Regulation 1935/2004 • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Nominal Volume: | 500ml |
| Brimful Volume: | 540ml |
| Method of Construction: | Single piece, injection moulded |
| Height: | 106mm |
| Lid Diameter: | 94mm |
| Nominal Weight: | 83g |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | N/A |
| Number of Bottles per Box: | 37 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |

ATEX:

This product has not been ATEX tested.

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-May-19