

Product Specification

Specification Number

8404H-1000S

Revision

02

Part No.

8404H-1000S

Description

Wide Mouth HDPE Bottle 1000ml - Sterile

Materials of Construction:

Bottle	HDPE	(conforms to FDA 177.1520 & USP Class VI)
Cap	PP	(conforms to FDA 177.1520 & USP Class VI)



Manufacturing Environment:

Controlled environment within factory

Packing Environment:

Class 100,000 Cleanroom

Nominal Volume:

1000ml

Brimful Volume:

1075ml

Nominal Height:

201mm

Diameter:

92mm

Internal Neck Diameter:

52mm

Weight:

100g

Individually Bagged:

Yes (in a PE/PET laminate film bag)

Pack Size:

30 bottle per box

Method of Sterilisation:

Gamma Irradiation (25-45 kGy)

BSE/TSE Free:

Yes

Supplied with lid screwed on:

Yes

Specification Approved By:

Signature:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-May-19