

Shakers, mixers and stirrers » MMS-3000 mini magnetic stirrer

MMS-3000 mini magnetic stirrer


Compact mini magnetic stirrer for routine laboratory procedures.

- Adjustable speed control: 0 to 3000rpm
- Maximum stirring volume: 20 litres



- Applications:**
- PH metering, extraction and dialysing with small quantities of substances
 - Any requirement for stirring solutions up to 20 litres

MMS-3000 mini magnetic stirrer – Technical specifications

		MMS-3000
		Mini magnetic stirrer
		 <p>h: 75mm d: 230mm w: 185mm weight: 1.5kg</p>
Operating temperature range	°C	+4 to 40
Speed	rpm	0 to 3000
Liquid stirring viscosity	mPa	Up to 1170
Stirring volume	L	20
Working plate size	mm	ø 160
External power supply		Input AC 120-230V, 50/60Hz Output DC 12V
Input voltage	V dc	12
Power consumption	W	3 (0.3A)