


## BagLight® PolySilk®

Blender bag

Ref 131 025 (100 mL) - Ref 132 025 (400 mL) - Ref 132 200 (2000 mL)  
Ref 133 025 (3500 mL)



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagLight® PolySilk® is a blender bag for all microbiological uses in labs.

### TECHNICAL SPECIFICATIONS

- Non-filter bag
- PolySilk®, polyolefin complex
- Flexible and transparent
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenization of the sample
- Compatible with any lab blender
- Approved for food contact: Regulation (EC) No 1935/2004
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Average shelf-life: 25 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 100/400/2000/3500 mL
- In compliance with ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

### BagLight® PolySilk® 100

Ref 131 025

- Max blending volume: 100 mL
- Optimal blending volume: 5 - 50 mL
- Bag dimensions: 110 x 200 mm
- Box dimensions: 35 x 24 x 9 cm, weight: 1.7 kg
- Pack of 25
- Box of 500

### BagLight® PolySilk® 400

Ref 132 025

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4 kg
- Pack of 25
- Box of 500

## BagLight® PolySilk® 2000

Ref 132 200

- Max blending volume: 2000 mL
- Optimal blending volume: 400 - 1500 mL
- Bag dimensions: 250 x 380 mm
- Box dimensions: 47.5 x 28.5 x 13.5 cm, weight: 9.6 kg
- Pack of 25
- Box of 500

## BagLight® PolySilk® 3500

Ref 133 025

- Max blending volume: 3500 mL
- Optimal blending volume: 400 - 3500 mL
- Bag dimensions: 380 x 560 mm
- Box dimensions: 44.5 x 25 x 35 cm, weight: 16 kg
- Pack of 25
- Box of 250