

Multi-Channel Solution Reservoirs

- Graduated inner walls
- Micro-trough™ U-shaped bottoms
- Pouring spouts on all 4 corners
- Sterile and disposable



	Cat. No.	Description	Qty
25mL	P8025	Reservoir, 25mL, bulk pack	100
	P8025-1S	Reservoir, 25mL, sterile, individually wrapped	100
	P8025-5S	Reservoir, 25mL, sterile, 5/sterile bag	200
50mL	P8050	Reservoir, 50mL, bulk pack	100
	P8050-1S	Reservoir, 50mL, sterile, individually wrapped	100
	P8050-5S	Reservoir, 50mL, sterile, 5/sterile bag	200
50mL Div	P8051	Reservoir, 50mL, 8 PLUS 4™ divided, bulk pack	100
	P8051-1S	Reservoir, 50mL, 8 PLUS 4™ divided, sterile, individually wrapped	100
	P8051-5S	Reservoir, 50mL, 8 PLUS 4™ divided, sterile, 5/sterile bag	200
100mL	P8100	Reservoir, 100mL, bulk pack	100
	P8100-1S	Reservoir, 100mL, sterile, individually wrapped	100
	P8100-5S	Reservoir, 100mL, sterile, 5/sterile bag	200

Standard 4-channel, 8-channel, 12-channel and various special purpose pipettes are all compatible with these disposable, polystyrene reservoirs. To facilitate easy, dripless pour-off of unused reagents and solutions, all versions have indented "spouts" in each corner. Each size is formed with a MicroTrough™ U-shaped bottom to minimize waste and reagent consumption. Packaged in notched, easy-open, tear strip bags.



ASPIR-8™ Reservoirs

- Specially designed for 8-Channel Pipettes
- 75mm Inner Chamber
- Conserve valuable reagents



Our new ASPIR-8™ is the ideal reservoir for any lab utilizing 8-channel pipettes. The strict inner dimensions are designed specifically for 8-channels, thus conserving space and reagent compared to its longer 12-channel counterparts that require 33% more reagent. This narrower design minimizes reagent exposure and usage. ASPIR-8 is the only reservoir that is tailor-made for your 8-channel instrument. Packaged in notched, easy-open, tear strip bags.

	Cat. No.	Description	Qty
5mL	P7005	ASPIR-8™ Reservoir for 8-Channel Pipettes, 5mL, bulk pack	300
	P7005-1S	ASPIR-8™ Reservoir for 8-Channel Pipettes, 5mL, sterile, individually wrapped	100
	P7005-5S	ASPIR-8™ Reservoir for 8-Channel Pipettes, 5mL, sterile, 5/sterile bag	200
10mL	P8010	ASPIR-8™ Reservoir for 8-Channel Pipettes, 10mL, bulk pack	300
	P8010-1S	ASPIR-8™ Reservoir for 8-Channel Pipettes, 10mL, sterile, individually wrapped	100
	P8010-5S	ASPIR-8™ Reservoir for 8-Channel Pipettes, 10mL, sterile, 5/sterile bag	200
25mL	P7025	ASPIR-8™ Reservoir for 8-Channel Pipettes, 25mL, bulk pack	100
	P7025-1S	ASPIR-8™ Reservoir for 8-Channel Pipettes, 25mL, sterile, individually wrapped	100
	P7025-5S	ASPIR-8™ Reservoir for 8-Channel Pipettes, 25mL, sterile, 5/sterile bag	200