

Human Dental Pulp Mesenchymal Stem Cells, Cryopreserved, 1.10^6 cells



Product sheet, catalog n° CTICC1.8.1 – Ref 169371

General Information

- **Organism:** Human (*Homo sapiens*)
- **Tissue:** Dental Pulp
- **Cell Type:** Dental Pulp Mesenchymal Stem Cells from single donor
- **Gender:** Male or Female (see Certificate of Analysis)
- **Age:** 18 to 99 years old (see Certificate of Analysis)
- **Donor type:** (see Certificate of Analysis)



Conservation:
 $T \leq -150^{\circ}\text{C}$



Biosafety level 1

Cell Characteristics

- **Cell properties:** Adherent
- **Morphology:** Bipolar, spindle-shaped, fibroblast-like morphology
- **Isolation:** Enzymatic dissociation
- **Cell viability:** Minimum 70% viability when thawed from cryopreservation
- **Cell conditioning:** Supplied as vials of 1.10^6 cells
- **Cryopreservation medium:** Frozen with 90% serum-free cryopreservation medium + 10% DMSO
- **Storage condition:** Liquid nitrogen
- **Batch specific information:** Included in the Certificate of Analysis

Safety and Quality Control

- **Biosafety level:** 1
- **Contamination:** Use mandatory laboratory protection and handle with care tissues and cells derived from human samples to avoid any contamination of the operator
- **Viral testing:** Negative for HIV, HBV, HCV
- **Sterility testing :** Negative for mycoplasma, bacteria and yeasts

Handling upon delivery and storage

- Check that the containers are intact and free of damage
- If not used immediately, place the vials at -150°C or below upon delivery

Thawing

1. Add 13 ml of PBS solution to a 15 ml conical tube, and warm in a water bath to 37°C
2. Thaw cryovial by swirling in a water bath at 37°C . As soon as thawed, start next step
3. Add the cell suspension to the warmed PBS in sterile conditions
4. Spin the tube at 250 g for 7 minutes to pellet the cells
5. Resuspend the cells in the appropriate volume of medium before using them

Associated products

- CTIOS.1.1: Adult Human Tooth, Fresh
- CTICC1.8.2: Human Oral Mucosa Fibroblasts, Cryopreserved, 1.10^6 cells
- CTICC1.8.3: Human Oral Mucosa Epithelial Cells, Cryopreserved, 1.10^6 cells

Provisions

- Cells and tissues are intended for **research use only** and shall not be used for human trials, animal trials, or diagnostics
- **Consent:** the original tissues have been obtained after informed consent of the patient under the provisions required by French Law
- **Primary Human cells** are not immortalised cell lines and may not be continually subcultured