


Human Placenta, Full Thickness, Flash Frozen, 10 cm³



Product sheet, catalog n° CTINT.1.4 – Ref : 169388

General Information

- **Organism:** Human (*Homo sapiens*)
- **Tissue:** Placenta
- **Location:** Whole Placenta (see Certificate of Analysis)
- **Gender:** Male or Female (see Certificate of Analysis)
- **Age:** New born (see Certificate of Analysis)
- **Donor type:** (see Certificate of Analysis)



Tissue Characteristics

- **Size:** 10cm³
- **Type of preparation:** Flash Frozen, unprocessed sample
- **Specificities:** For specific requirements as pathology, etc. (see Certificate of Analysis)

Safety and Quality Control

- **Biosafety level:** 1
- **Contamination:** Use mandatory laboratory protection and handle with care tissues and cells derived from human samples to avoid any contamination of the operator.
- **Viral testing:** Negative for HIV, HBV, HCV

Handling upon delivery and storage

- Check that the containers are intact and free of damage
- If not used immediately, place the vials at -80°C or below upon delivery

Recommended applications

- **Type of analysis:** Histological, Immunofluorescence, Molecular biology
- **Type of testing:** Active components, Toxicology, Biomarkers studies, Dynamic and Mechanical properties
- **Fields of study:** Anti-inflammatory activity, Drug development, Pregnancy Disorders, Healing and Repair, Immunogenicity

Associated products

- **CTINT.1.8:** Human Amniotic membrane, Flash Frozen, 10cm²
- **CTINT.1.12:** Human Umbilical Cord, Flash Frozen, 10cm
- **CTICC.1.12:** Human Amniotic Membrane Mesenchymal Stem Cells, Cryopreserved, 1.10⁶ cells

Provisions

- Cells and tissues are intended for **research use only** and shall not be used for human trials, animal trials, or diagnostics.
- **Consent:** the original tissues have been obtained after informed consent of the patient under the provisions required by French Law.