

Human Peripheral Blood Plasma, Fresh or Frozen, 2mL



Product sheet, catalog n° CTIPB.1.15 – Ref : 169421

General Information

- **Organism:** Human (*Homo sapiens*)
- **Tissue:** Peripheral Blood
- **Gender:** Male or Female (see Certificate of Analysis)
- **Age:** 18 to 99 years old (see Certificate of Analysis)
- **Donor type:** (see Certificate of Analysis)



Tissue Characteristics

- **Volume:** 2mL
- **Type of preparation:** Fresh or frozen, unprocessed sample
- **Specificities:** For specific requirements as pathology, etc. (see Certificate of Analysis)

Safety and Quality Control

- **Biosafety level:** 1
- **Contamination:** Use mandatory laboratory protection and handle with care tissues and cells derived from human samples to avoid any contamination of the operator.
- **Viral testing:** Negative for HIV, HBV, HCV

Handling upon delivery and storage

- Check that the containers are intact and free of damage
- If not used immediately, place the vials at +4°C (if fresh) or -80°C (if frozen) upon delivery.

Recommended applications

- **Type of analysis:** Molecular biology, Immunology, Haematology
- **Type of testing:** Active components, ADME-Toxicology, Biomarkers studies, Pharmaceuticals
- **Fields of study:** Anti-inflammatory activity, Drug development, Immunogenicity, Cell isolation, Cell Therapy, Regenerative Medicine, Transplantation

Associated products

- CTIPB.1.8: Human Peripheral Blood Serum, 2mL
- CTIPB.1.16: Human Peripheral Blood Plasma, 10mL
- CTIPB.1.18: Human Peripheral Blood Buffy Coat, 2mL

Provisions

- Cells and tissues are intended for **research use only** and shall not be used for human trials, animal trials, or diagnostics.
- **Consent:** the original tissues have been obtained after informed consent of the patient under the provisions required by French Law.
- **Primary Human cells** are not immortalised cell lines and may not be continually subcultured.