

# Adult Human Peripheral Blood Monocytes, CD14+, Cryopreserved, 1.10<sup>7</sup> cells



Product sheet, catalog n° CTICC.IMPB.8 – Ref : 169434

## General Information

- **Organism:** Human (*Homo sapiens*)
- **Tissue:** Peripheral Blood
- **Cell Type:** Peripheral blood mononuclear cell
- **Gender:** Male or Female (see Certificate of Analysis)
- **Age:** 18 to 99 years old (see Certificate of Analysis)
- **Donor type:** (see Certificate of Analysis)



Conservation:  
T ≤ -150 °C



Biosafety level 1

## Cell Characteristics

- **Cell properties:** Suspension
- **Morphology:** Spherical
- **Cell surface marker:** CD14+
- **Isolation:** Density gradient separation Ficoll Plaque plus + Immunomagnetic selection
- **Cell viability:** Minimum 80% viability when thawed from cryopreservation
- **Cell conditioning:** Supplied as vials of 1.10<sup>7</sup> cells
- **Cryopreservation medium:** Frozen with 90% serum-free cryopreservation medium + 10% DMSO
- **Storage condition:** Liquid nitrogen
- **Batch specific information:** Included in the Certificate of Analysis

## Safety and Quality Control

- **Biosafety level:** 1
- **Contamination:** Use mandatory laboratory protection and handle with care tissues and cells derived from human samples to avoid any contamination of the operator.
- **Viral testing:** Negative for HIV, HBV, HCV

## Handling upon delivery and storage

- Check that the containers are intact and free of damage
- If not used immediately, place the vials at -150°C or below upon delivery

## Thawing

1. Add 13 ml of PBS solution to a 15 ml conical tube, and warm in a water bath to 37 °C
2. Thaw cryovial by swirling in a water bath at 37 °C. As soon as thawed, start next step
3. Add the cell suspension to the warmed PBS in sterile conditions
4. Spin the tube at 250 g for 7 minutes to pellet the cells
5. Resuspend the cells in the appropriate volume of medium before using them

## Associated products

- **CTICC.IMPB.1:** Adult Human Peripheral Blood PBMC, Cryopreserved, 10.10<sup>6</sup> cells
- **CTICC.IMPB.3:** Adult Human Peripheral Blood B-Cells, CD19+, Cryopreserved, 1.10<sup>6</sup> cells
- **CTICC.IMPB.5:** Adult Human Peripheral Blood T-Cells, CD3+, Cryopreserved, 1.10<sup>7</sup> cells

## Provisions

---

- Cells and tissues are intended for **research use only** and shall not be used for human trials, animal trials, or diagnostics.
- **Consent:** the original tissues have been obtained after informed consent of the patient under the provisions required by French Law.
- **Primary Human cells** are not immortalised cell lines and may not be continually subcultured.

DOMINIQUE DUTSCHER SAS