

# 96-Well Semi-Skirted Plate, *Newly Improved*



➤ Extremely well-manufactured plate *is flat, and stays flat*

➤ Manufactured to SBS standards; fits a variety of thermal cyclers and sequencers

➤ Ultra thin-well walls for maximum heat transfer

➤ Holes between wells minimize vacuum effect in cycling

➤ Black printed alpha-numeric provide easy sample identification

➤ 100% quality control

➤ Certified *free* of RNase/DNase, human DNA, and pyrogens

➤ Compatible with Sorenson plate seals and cap strips

➤ Bar coding available

