

TECHNICAL DATA SHEET

Part Number	17836-50G-HON
Description	Ammonium acetate

Technical specifications :

- ✓ ammonium acetate for hplc =99,0% 50g
- ✓ CAS N°631-61-8
- ✓ Chemical purity ≥99%
- ✓ Molecular weight 77.08 g/mol
- ✓ Molecular formula C₂H₄O₂.H₃N
- ✓ Density 1.170 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form solid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 50g

- ✓ P Statements P280

Sales unit :

50g Ammonium acetate
