

TECHNICAL DATA SHEET

Part Number	17843-50G-HON
Description	Ammonium formate

Technical specifications :

- ✓ ammonium formate, for hplc =99,0% 50g
- ✓ CAS N°540-69-2
- ✓ Chemical purity ≥99%
- ✓ Molecular weight 63.06 g/mol
- ✓ Molecular formula CH₅NO₂
- ✓ Density 1.266 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form crystalline

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 50g

- ✓ GHS pictograms GHS07
- ✓ H Statements H319
- ✓ P Statements P280

Sales unit :

50g Ammonium formate
