

TECHNICAL DATA SHEET

Part Number	17844-250G-HON
Description	Sodium phosphate monobasic

Technical specifications :

- ✓ sodium phosphate monobasic, for hplc =99,0% 250g
- ✓ CAS N°7558-80-7
- ✓ Chemical purity ≥99%
- ✓ Molecular weight 119.98 g/mol
- ✓ Molecular formula H₂NaO₄P
- ✓ Density 2.36 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form crystalline

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 250g

- ✓ P Statements P280

Sales unit :

250g Sodium phosphate monobasic
