



Brand: Honeywell
Product: 179973

Acetone Laboratory Reagent, 99.5%

Parameter	Specification	Units
IR spectrum	complying	
assay (GC)	min. 99.50	%
refractive index (n 20/D)	1.357-1.361	
non-volatile matter	max. 0.002	%
water (Karl Fischer)	max. 0.5	%

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.