



ClearLine®

ClearLine® Adhesive Lab Tape

Rolls of coloured adhesive tape made for use in laboratories

Features :

- ✓ Robust, resistant to water, oils and acids
- ✓ Can be used in a freezer at -80°C and autoclaved 15min at +121°C
- ✓ Can be written on with a pen, pencil or permanent marker
- ✓ Adheres to most surfaces : glass, plastic and metals
- ✓ Removes and cleans easily and quickly

- ✓ Available in 12m and 55m -long rolls , and with a 13, 19 or 25mm width

TECHNICAL DATASHEET

REFERENCE :

Rainbow Packs

	Content (in rolls)							Cat No
	White	Yellow	Red	Orange	Green	Blue	Pink	
24 rolls 12m x 13mm	5	4	3	3	3	3	3	188007
16 rolls 12m x 19mm	4	2	2	2	2	2	2	188015
12 rolls 12m x 25mm	2	2	2	1	2	1	2	188023

Adhesive Lab Tape – 12m long

Tape width	13 mm		19 mm		25 mm	
	6 rolls	1 roll	4 rolls	1 roll	3 roll	1 roll
White	188000	188000A	188008	188008A	188016	188016A
Yellow	188001	188001A	188009	188009A	188017	188017A
Green	188002	188002A	188010	188010A	188018	188018A
Red	188003	188003A	188011	188011A	188019	188019A
Orange	188004	188004A	188012	188012A	188020	188020A
Blue	188005	188005A	188013	188013A	188021	188021A
Pink	188006	188006A	188014	188014A	188022	188022A
Turquoise	188048	188048A	188050	188050A	188052	188052A
Fuchsia	188049	188049A	188051	188051A	188053	188053A

Adhesive Lab Tape – 55m long

Tape width	13 mm		19 mm		25 mm	
	6 rolls	1 roll	4 rolls	1 roll	3 rolls	1 roll
White	188024	188024A	188031	188031A	188038	188038A
Yellow	188025	188025A	188032	188032A	188039	188039A
Green	188026	188026A	188033	188033A	188040	188040A
Red	188027	188027A	188034	188034A	188041	188041A
Orange	188028	188028A	188035	188035A	188042	188042A
Blue	188029	188029A	188036	188036A	188043	188043A
Pink	188030	188030A	188037	188037A	188044	188044A
Turquoise	188054	188054A	188056	188056A	188058	188058A
Fuchsia	188055	188055A	188057	188057A	188059	188059A

Tape Dispensers

188045	ClearLine Lab tape dispenser – for 55m rolls
188046	ClearLine Lab tape dispenser – for 12m rolls
188047	Single-roll ClearLine Lab tape dispenser – for 12m rolls