


TECHNICAL DATA SHEET

Part Number	191018
Description	Centrifuge Tubes 50 ml



Benefits :

- ✓ Tube 50 ml, conical - max. g-Force : 15'500
- ✓ Manufactured with highest grade of raw material fulfilling requirements of USP Class VI
- ✓ white writing circle, affixed on the screw cap, for easy labeling with any writing instrument
- ✓ caps that cannot be over-tightened
- ✓ gas and aerosol tight cap
- ✓ graduations all the way down to 0.1 ml / 0.5 ml respectively
- ✓ convenient racks available that are autoclavable and are collapsible for storage when not in use

Technical specifications :

- ✓ Rack "L", 20 x 50 ml
- ✓ Material : Polypropylene

Sales unit :

200 Tubes