

TECHNICAL DATA SHEET

KSS-67810-LT

PRODUCT CODE

KSS-67810-LT

PACKAGING

Cases of 100 containing 5 packs of 20 units

SHELF LIFE

15 months (upon arrival store between 2°C & 8°C)

PRODUCT DESCRIPTION

SANI-SPONGE kits are designed to collect environmental samples and used to detect microbiological contaminants such as Listeria, Salmonella, E. Coli and other food-borne pathogens on food preparation surfaces. These kits are generally used in the food, environmental and cosmetic industries. The SANI-SPONGE kit can be sold with a dry sponge or a sponge soaked in various solutions for transportation.

PRODUCT SPECIFICATIONS

Closure type	2 round wires closure system with safety tabs
Closure material	Steel wires covered with FDA approved pressure sensitive tape with a backing of polypropylene or PVC.
Printed	Sani Sponge printed
MSDS / Formulation	Available for the solution only, upon request
Sterilisation	Sterility of product is ensured by gamma irradiation. Validation tests based on ISO 11137-1 and 11137-2 are performed in our facility to guaranty sterility of our processes.
Dimensions & volume	7.5" x 11.5" (191 mm x 292 mm) capacity 1330 mL (45 oz)
Thickness of material	3.5 mil (89 microns)
Solution	Lethen broth
Volume of solution	10 mL
Composition of sponge	Polyurethane
FDA regulation	All TWIRL'EM sterile sampling bags comply to the applicable regulations of the Food and Drug Administration (21 CFR 177.1520 (c)3.1a, (c)3.2a and 178.2010.
EU regulation	All TWIRL'EM sterile sampling bags comply where applicable to EU No 10/2011 regulation for use in contact with food, subject to certain types of migration limitations and labeling requirements.
CFIA	All TWIRL'EM sterile sampling bags are also a solution which could be used in the CFIA Prevention Control Plan inspired by HACCP principles, a Canadian federal initiative from the Safe Food for Canadians Regulations (SFCR).

T : 450-649-7343

F : 450-649-3113

www.labplas.com
 service@labplas.com
 Update : 04/11/2021