

FILTRA-BAG

Once at the laboratory, solid food samples must be separated and homogenized. Labplas FILTRA-BAG is perfect for use in the homogenizer.

Labplas FILTRA-BAG blender bags are designed to simplify taking an aliquot when working with samples which contain large amounts of residue and/or semi-solid/solid substances. When lacing an aliquot onto growth media, it's very important that the volume of the aliquot be accurate and that it be free of solid particles. This prevents pipette blockage as well as uneven plate distribution.

Product Advantages:

- Made with a polyethylene/Pet outer shell which practically eliminates the possibility of piercing the bag wall during blending.
- The dividing filter membrane has 1840 holes per square inch with a pore size of approximately 330 microns. This pore size has proven to be very effective for the majority of applications and allows for optimal solution and bacterial flow between the FILTRA-BAG® compartments during blending.



FILTRA-BAG with Open Top - Clear

CODE	MIL. IN	MICR.	IN X IN	MM X MM	OZ	ML	PACKAGING
SCTO7012A - For 400 series blenders	4.0	101	7 x 12	178 x 305	55	1650	Cs/400 (4x100)

FILTRA BAG with Tear-Off Protection Strip - Clear

CODE	MIL. IN	MICR.	IN X IN	MM X MM	OZ	ML	PACKAGING
SCT-6090A - For 400 series blenders	4.0	101	6 x 9	152 x 229	24	710	Cs/400 (4x100)
SCT-7012A - For 400 series blenders	4.0	101	7 x 12	178 x 305	55	1650	Cs/400 (4x100)
SCT-1015A - For 3500 series blenders	4.0	101	10 x 15	254 x 381	92	2720	Cs/400 (4x100)
SCT-1515A - For 3500 series blenders	4.0	101	15 x 15	381 x 381	138	4000	Cs/200 (2x100)
SCT-1520A - For 3500 series blenders	4.0	101	15 x 20	381 x 508	405	12L	Cs/200 (2x100)

FILTRA BAG with - Safety Tabs - Clear - Closure with 1 round wire & 1 flat wire

CODE	MIL. IN	MICR.	IN X IN	MM X MM	OZ	ML	PACKAGING
EFT-6090A - For 400 series blenders	4.0	101	6 x 9	152 x 229	24	710	Cs/400 (4x100)
EFT-7012A - For 400 series blenders	4.0	101	7 x 12	178 x 305	55	1650	Cs/400 (4x100)
EFT-1015A - For 3500 series blenders	4.0	101	10 x 15	254 x 381	92	2720	Cs/400 (4x100)
EFT-1515A - For 3500 series blenders	4.0	101	15 x 15	381 x 381	138	4000	Cs/200 (2x100)
EFT-1520A - For 3500 series blenders	4.0	101	15 x 20	381 x 508	405	12L	Cs/200 (2x100)