

# Product Specification — CO-RE Tip

## P/N and Product description

---

- P/N 194050
- 5ml 1T Rack CO-RE Tip, non-sterile, black, for Microlab STAR Line



## Material

---

- Polypropylen
- conductive

## Packaging & processing

---


- 4 tips/tray, 24 trays/case, 96 tips/case
- Color-code label on tip tray: green
- Barcode label contains: unique identification number, Lot No., production date, cavity no.
- Barcode type: Code 128 subset B or C
- Tips usable on CO-RE single channel and 5000µl pipetting channels



## Certification / Purity grade

---

- ISO Class 8 Clean Room production, compliant with ISO 9001, GMP and FDA
- CE IVD conform



- DNase-/RNase-free
- Pyrogen-free (endotoxin-free)
- PCR inhibitor-free

## Labeling

---

Each box of disposable tips provides a label which indicates a barcode, part number (REF) and Lot No.

## Measurement

---

- Tip length: 118mm
- Inner diameter tip orifice: 1.05mm

**HAMILTON** 

Web: [www.hamiltonrobotics.com](http://www.hamiltonrobotics.com)

Email: [infoservice@hamiltonrobotics.com](mailto:infoservice@hamiltonrobotics.com)

United Kingdom & Ireland

Tel: +44 (0)121-717-0199

China

Tel: +86-21-6164-6567

France

Tel: +33 (01) 69751616

Italy

Tel: +39-39-689-33-93

Denmark, Norway,  
Sweden, Finland

Tel: +46 (0) 8 410 273 73

Germany, Switzerland,  
Austria, Benelux

Tel: +49 (089) 552649-0